
Venus Documentation

Release 0.1

TChapman

Feb 23, 2024

CONTENTS

1	Venus	3
1.1	General Steps	3
1.2	ML_STAR	4
1.3	Custom Dialog Steps	5
2	Error Codes & Debugging	7
2.1	Major IDs that have been identified	7
2.2	Minor IDs that have been identified	11
2.3	Internal Error Codes which have been identified (presumably is the same as SpecificErrorId?)	14
3	ASWGlobal	19
4	Aliquot for Easy Steps v2	21
5	Aliquot with Easy Steps	23
6	Array (from HSLExtensions)	25
7	ArrayTools	31
8	ColourSchemeLibrary	35
9	Convert File To ASCII	37
10	DebugAssist	39
11	EditFileAttributes	41
12	Error Report Library	43
13	Error Simulator	45
14	Framework (from HSLExtensions)	51
15	HSLDeckVisualize	53
16	HSLFillLib	55
17	HSLFileLibEx	61
18	HSLStatistics	63
19	HSLStrLib	67

20	HSLZipLib	77
21	HSL_SeqDailyTools	79
22	If_And_If_Or	81
23	Labware Properties	83
24	Lookup	89
25	Pipetting (from HSLExtensions)	91
26	STAR_Channel_Tools	95
27	StrTokenize	113
28	String (from HSLExtensions)	115
29	Tools Library	119
30	TraceLevel	121
31	Windows (from HSLExtensions)	129
32	Zero uL Scanner	131
	Index	133

Venus is a coding language wrapper used in the coding of Hamilton Microlab STARS and similar devices. This is the consolidated unofficial documentation for Venus 4 as well as its libraries, compiled from official sources as well as personal troubleshooting.

This is not an official Hamilton resource

Note: This project is under active development. If any issues or mistakes are found, please comment or contact me directly at tarunchapman@hotmail.com. General place to download libraries: <https://github.com/theonetrueerd/VenusPackages>

VENUS

This section aims to cover the *functions* that Venus can run. Documentation on more general principles of Venus, including things like the liquid editor and labware editor, will be found separately.

NB: This section is probably coming last as out of all of the Venus resources it has the best documentation.

The Venus functions are split into a few main categories:

- General Steps
- ML_STAR
- Custom Dialog Steps

1.1 General Steps

The general steps are steps which are not instrument-specific (and thus not hardware-related), and aren't involved in custom dialog creation. They are the following:

- `Comment()`
- `Assignment()`
- `Assignment with Calculation()`
- `Loop()`
- `Loop: Break()`
- `If, Else()`
- `Array: Declare / Set Size()`
- `Array: Set At()`
- `Array: Get At()`
- `Array: Get Size()`
- `Array: Copy()`
- `Sequence: Get Current Position()`
- `Sequence: Set Current Position()`
- `Sequence: Get End Position()`
- `Sequence: Set End Position()`
- `Adjust Sequences()`

- File: Open()
- File: Read()
- File: Write()
- File: Set Position()
- File: Close()
- Timer: Start()
- Timer: Wait For()
- Timer: Read Elapsed Time()
- Timer: Restart()
- User Input()
- User Output()
- Shell()
- Set Event()
- Wait for Event()
- Return()
- Abort()
- Error Handling by the User()
- Begin Parallel()
- End Parallel()

1.2 ML_STAR

The ML_STAR steps are steps which relate to the ML_STAR and interact with its hardware specifically, doing things such as pipetting, moving channels, etc. They are split into two categories; the four Smart Steps and the remaining Single Steps. They are as follows:

- 1000uL Channel Aspirate()
- 1000uL Channel Dispense()
- iSWAP Transport()
- 1000uL Channel CO-RE Grip Transport()
- 1000uL Channel Tip Pick Up (Single Step)()
- 1000uL Channel Aspirate (Single Step)()
- 1000uL Channel Dispense (Single Step)()
- 1000uL Channel Dispense on the Fly (Single Step)()
- 1000uL Channel Tip Eject (Single Step)()
- 1000uL Channel Get Last Liquid Level (Single Step)()
- 1000uL Channel Aspirate 2nd Phase (Single Step)()
- Initialize (Single Step)()

- Lock/Unlock Front Cover (Single Step)()
- iSWAP Get Plate (Single Step)()
- iSWAP Place Plate (Single Step)()
- iSWAP Move Plate (Single Step)()
- iSWAP Open Gripper (Single Step)()
- iSWAP Close Gripper (Single Step)()
- iSWAP Get First Plate Position (Single Step)()
- iSWAP Park (Single Step)()
- 1000uL Channel CO-RE Grip Get Plate (Single Step)()
- 1000uL Channel CO-RE Grip Place Plate (Single Step)()
- 1000uL Channel CO-RE Grip Move Plate (Single Step)()
- 1000uL Channel Move To Position (Single Step)()
- Wait for TADM Upload (Single Step)()

1.3 Custom Dialog Steps

The custom dialog steps only has a single step called custom dialog, which helps in the creation of more personalised versions of the standard User Input and User Output dialogs.

- Custom Dialog()

ERROR CODES & DEBUGGING

I will try compile as many different error codes and their explanations here, as well as potential causes and solutions to deal with them as and when they pop up. If there are any you wish to add, please message me on my discord “theonetrue nerd” or email me at “tarunchapman@hotmail.com”. Notably, the error codes that show up in the trace look slightly different to the ones presented here. If the trace error code looks like this: (0x12 - 0x3 - 0x45) [which corresponds to (MinorID - MajorID - SpecificErrorID)], then the translated error code will look like this: 0xa3120045. The error code should always be 10 digits long, with zeroes padding the bit between the 2 and 45.

The Major ID is usually assigned in the library which is causing the error; for example in HSLUtilLib2 there is a function called “Error” in which there is a static const variable called MajorID. The minor ID corresponds to the general type of error; in HSLUtilLib2 “0x01” is used to refer to a general runtime error.

At least I think? I’ve found a bit of inconsistency so will keep digging.

- Runtime Errors
- Syntax Errors
- Unsorted Errors

2.1 Major IDs that have been identified

NB: The codes are in hexadecimal so (for example) 31 = 1F

- 0: None
- 2: HxCfgFil
- 3: HxTrcFil
- 4: HxParams
- 5: HxLabwr
- 6: HxPlumb
- 7: HxM4Comm
- 8: HxM4Inst
- 9: HxM4PPtr
- 10: HxPm4Met
- 11: HxPm4Cfg
- 12: HxPm4Run
- 13: HxPhnx

- 14: HxMovPrb
- 15: HxSsInst
- 16: HxDeckEd
- 17: HxRun
- 18: HxReg
- 19: HxM4MetsDat
- 20: HxM4MetsCfg
- 21: HxP96P4Wz
- 22: HxP96M4MetsRun
- 23: HxM4PPtrCfg
- 25: HxM4LoadCfg
- 26: HxM4PrbInLabwr
- 27: HxCfgEd
- 28: HxPm4Wz
- 29: HxWrkFil
- 30: HxM4PPtrCfg2
- 31: HxM4CComm
- 32: HxM4PPtrPars
- 33: HxM4Transfer
- 34: HxHslParser
- 35: HxHslExecutor
- 36: HxHslRctl
- 37: HxFdxProtocol
- 38: HxSerial
- 39: HxTextMet
- 40: HxGruCommand
- 41: HxGruInst
- 42: HxCommand
- 43: HxSwapCommand
- 44: HxProSim
- 45: HxCompCommand
- 46: HxMetEdCompCmd
- 47: HxPrbInLw
- 48: HxMetEd
- 49: HxCommunication
- 50: HxProtocol

- 51: HxUsbComm
- 52: HxGruCompCmd
- 53: HxSampleTracker
- 54: HxGruLiquid
- 55: HxSecurity
- 56: HxGruLiquidEditor
- 57: HxStarMaintAndVer
- 58: HxSoftMaxPro
- 59: HxCytomat
- 60: HxM4Command
- 61: HxCarousel
- 62: HxMetChaining
- 63: HxAuditTrail
- 64: HxTrace
- 65: HxTraceView
- 66: HxView
- 67: HxWatchView
- 68: HxPowerWaveHt
- 69: HxInstrumentData
- 70: HxMethodCopy
- 71: HxBtiElx405AutoWasher
- 72: HxSecurityCom
- 73: HxStarBvsCommand
- 74: HxStarBvsConfig
- 75: HxReportConfig
- 76: HxServices
- 77: HxScheduleView
- 78: HxStCompCmd
- 79: HxStarData
- 80: HxUserManager
- 81: HxSchedCompCmd
- 82: HxFan
- 83: HxSysDeck
- 84: HxVSpin
- 85: HxElementCounter
- 86: HxStarConfig

- 87: HxConfigEditor
- 88: HxStarDevices
- 89: HxVSpinAccess2
- 90: HxEditSequenceDlg
- 91: HxStarBiotechMaintMet
- 92: HxLabwrcat
- 93: HxLabwrCatComponents
- 94: HxStarDynDilLib
- 95: HxMosquito
- 96: HxDatabase
- 97: HxVectorDatabase
- 98: HxTranslationSupport
- 99: HxUtilLib
- 100: HxReagentDisp
- 101: HxSpeVacuum
- 102: HxM384
- 103: HxGen5
- 104: HxBigBearShaker
- 105: HxSys3DView
- 106: HxImpactCmd
- 107: HxXRPLiteMC
- 108: NimbusFourProbe
- 109: TrackGripper
- 110: HxCoreDevices
- 111: HxCoreLiquid
- 112: HxCoreLiquidEditor
- 113: HxXRPLiteConfigurator
- 114: HxTcpIpBdzComm
- 115: MiroIncubatorCmd
- 116: GuavaLinkCmd
- 117: PowerSocket
- 118: HxCustomDialog
- 119: GlasColMixerCmd
- 120: EntryExit
- 121: ForteOctetCmd
- 122: NexusXPeelCmd

- 123: PHERAstar_module
- 124: HXETRACKCMD

2.2 Minor IDs that have been identified

NB: Codes are once again in hexadecimal. I'm slightly unsure on some of these

HxVectorDB:

- 0: UnexpectedError
- 1: InterfaceNotInitialized
- 2: CannotInitializedInterface
- 3: CannotGetTableSchema
- 4: CannotGetElementCounter
- 5: ErrorInWorkerThread
- 6: CannotGetProperty
- 7: CannotSetProperty
- 8: CannotGetTranslationTable
- 9: CannotAddAdditionalData
- 10: CannotGetAdditionalData
- 11: CannotDeleteAdditionalData
- 16: Worklist_CannotAddJobs
- 17: Worklist_CannotGetJobs
- 18: Worklist_CannotGetJobAdditionalData
- 19: Worklist_CannotRemoveJob
- 20: Worklist_CannotGetJobState
- 21: Worklist_CannotAddJobAdditionalData
- 22: Worklist_CannotDeleteJobAdditionalData
- 23: Worklist_CannotGetJob
- 24: Worklist_CannotSetJobState
- 48: Tracking_CannotStartRun
- 49: Tracking_CannotPauseRun
- 50: Tracking_CannotResumeRun
- 51: Tracking_CannotEndRun
- 52: Tracking_CannotAbortRun
- 53: Tracking_CannotInterruptRun
- 54: Tracking_CannotCreateDeck
- 55: Tracking_CannotGetDeckID

- 56: Tracking_CannotGetAllLabwareOnDeck
- 57: Tracking_CannotGetLoadStateOfLabware
- 58: Tracking_CannotGetElementID
- 59: Tracking_CannotGetLabware
- 60: Tracking_CannotGetLabwareLoadingTime
- 61: Tracking_CannotAssignLabwareToJob
- 62: Tracking_CannotAssignLabwareToJobs
- 63: Tracking_CannotGetLabwareBarcode
- 64: Tracking_CannotGetLabwareVolume
- 65: Tracking_CannotGetLabwareState
- 66: Tracking_CannotUpdateTADMCurveIDForVolumeMove
- 67: Tracking_CannotGetLabwareLastSourceBarcode
- 68: Tracking_CannotGetLabwareSourceBarcodeList
- 69: Tracking_CannotGetLabwareInitialValues
- 70: Tracking_CannotGetRunID
- 71: Tracking_CannotGetRun
- 72: Tracking_CannotGetRunState
- 73: Tracking_CannotGetUserRunState
- 74: Tracking_CannotSetUserRunState
- 75: Tracking_CannotGetRunActions
- 76: Tracking_CannotGetRunAction
- 77: Tracking_CannotAddRunAction
- 78: Tracking_CannotGetElementIDs
- 79: Tracking_CannotGetLabwareHierarchy
- 80: Tracking_CannotTrackActionLoad
- 81: Tracking_CannotTrackActionUnload
- 82: Tracking_CannotTrackActionMoveVolume
- 83: Tracking_CannotTrackActionMove
- 84: Tracking_CannotTrackActionWash
- 85: Tracking_CannotTrackActionIncubate
- 86: Tracking_CannotTrackActionSetBarcode
- 87: Tracking_CannotTrackActionSetLabwareState
- 88: Tracking_CannotTrackActionCustomAction

HxDatabase:

- 00: Global
- 01: Singleton

- 02: ResourceManager
- 03: Error
- 05: IHxDbManagement
- 06: IHxDbConfiguration
- 07: Utilities
- 16: IHxDbCommand
- 17: IHxDbCommandCollection
- 18: IHxDbConnection
- 19: IHxDbCreateParameterCollection
- 20: IHxDbCreateProcedureCommand
- 21: IHxDbCreateTableCommand
- 22: IHxDbDataReader
- 23: IHxDbParameter
- 24: IHxDbParameterCollection
- 25: IHxDbTransaction

[Need sorting]

- 10: HxEmail
- 20: HxErrorEvent
- 21: SendAddress
- 22: SendFlag
- 23: ExecuteFlag
- 24: ExecuteName
- 26: startApplication
- 27: ExecuteArgument
- 40: HxGeneralSettings()
- 42: RequiredPasswordLength
- 43: SystemName
- 44: SimulationOn
- 45: GetSequenceRGB
- 46: AskForSequenceNameAfterDrop
- 47: GetSound
- 48: SetSound
- 49: GetTimeout
- 50: SetTimeout
- 60: HxInstallation
- 62: GetFeatureDescription

- 63: ActivateFeature
- 64: InstallFeature
- 65: GetFeatureExpiryISODate
- 66: GetFeatureStatusText
- 67: LegalizeInstallation
- 68: GetFeatureNameFromId
- 69: GetFeatureDescriptionById
- 70: GetFeatureExpiryISODateById
- 71: GetFeatureStatusTextById
- 72: UninstallFeature

2.3 Internal Error Codes which have been identified (presumably is the same as SpecificErrorId?)

NB: Once again, these are in Hexadecimal and four digits - error code 10 is 000A

HxVectorDatabase:

- 0: UnexpectedError
- 1: InterfaceAlreadyInitialized
- 2: CannotDetermineAdditionalValueType
- 3: UnknownEnumValue
- 4: UnsupportedEnumValue
- 5: UnsupportedValueFromDB
- 6: InconsistentDB
- 7: NotImplemented
- 8: CannotConnectToDb
- 9: CannotDetermineOS
- 10: ValueOutOfRange
- 11: ExternalDatabaseServerNotSupportedInStandardVersion
- 12: FunctionNotSupportedInStandardVersion
- 13: FunctionNotSupportedOnRemoteDatabaseServer
- 14: CannotFindConfigFile
- 15: CannotFindSqlScript
- 16: CannotConvertValueToByte
- 17: CannotConvertValueToShort
- 18: CannotConvertValueToInt
- 19: CannotConvertValueToLong

- 20: CannotConvertValueToBool
- 21: CannotConvertValueToDouble
- 22: CannotConvertValueToString
- 23: CannotConvertValueToDateTime
- 32: CannotConvertValueToHxVectorDbJobState
- 33: CannotConvertValueToHxVectorDbValueType
- 34: CannotConvertValueToHxVectorDbActionState
- 35: CannotConvertValueToHxPars
- 36: CannotConvertValueToHxVectorDbLabwareHandling
- 37: CannotConvertValueToHxVectorDbLabwareLevel
- 38: CannotConvertValueToHxVectorDbLabwareState
- 39: CannotConvertValueToHxVectorDbStepType
- 40: BadParameterSupplied
- 41: CannotConvertValueToHxVectorDbRunState
- 42: CannotSetValueInConfigFile
- 43: CannotGetValueInConfigFile
- 44: CannotFindStringInStringTable
- 45: CannotConvertValueToHxVectorDbActionType
- 46: CannotConvertValueToHxVectorDbLabwareUsageType
- 47: CannotConvertValueToHxVectorDbRunAction
- 48: CannotDetachDatabase
- 49: CannotConvertValueToHxVectorDbSortingAlgorithm
- 50: CannotRenameExistingDatabaseFiles
- 51: CannotConvertValueToHxVectorDbActionType
- 52: SqlScriptChecksumVerificationFailed
- 4097: CannotFindJob
- 4098: CannotDeterminePhoenixVersion
- 4099: CannotDetermineCurrentUsername
- 4100: RunAlreadyStarted
- 4101: RunNotRunning
- 4102: RunNotPaused
- 4103: MultipleRunsWithSameGUIDDetected
- 4104: CannotUpdateInternalRunID
- 4105: CannotDetermineInstrumentIDForConfiguration
- 4106: InstrumentDuplicatesInDatabase
- 4107: InstrumentConfigurationDuplicatesInDatabase

- 4108: DeckAlreadyExistsForInstrument
- 4109: LabwareAlreadyExists
- 4110: LabwareDoesNotExist
- 4111: OnlySingleLabwareAccessAllowed
- 4112: UnknownDeckID
- 4113: IllegalLabwareHandlingCombination
- 4114: CannotExtractLabwareName
- 4115: CannotMixActionsOfDifferentActionTypes
- 4116: IfNotExistsCreateNotAllowed
- 4117: IfExistsRemoveNotAllowed
- 4118: IfNotExists_ErrorNotAllowed
- 4119: CannotExtractBaseName
- 4120: CannotNormalizeLabwareAccessName
- 4121: CannotLinkLabware
- 4122: CannotExtractInstrumentName
- 4123: CannotSplitLabwareAccessName
- 4124: CannotDetermineLabwareStatePriority
- 4125: CannotConvertActionStateToLabwareState
- 4126: WrongNumberOfLabwareForAction
- 4127: CannotFindRun
- 4128: CannotFindAction
- 4129: CannotFinishPreWork
- 4130: CannotCreateAction
- 4131: CannotLinkLabware
- 4132: CannotAddDetail
- 4133: CannotAddAdditionalData
- 4134: CannotUpdateLabwareData
- 4135: IfNotExistsIgnoreActionNotAllowed
- 4136: DeckDoesNotExistForInstrument
- 4137: SequenceNotValidAtIndex
- 4138: CannotAddAdditionalData
- 4139: CannotGetAdditionalData
- 4140: CannotDeleteAdditionalData
- 4141: CannotFindRunAction
- 4142: CannotFindInstrument
- 4143: SequenceNotValid

- 4144: ConnectedContainerNotAllowed
- 4145: LabwareOrLabwareTypeNotValid
- 4146: LabwareTypeDoesNotExist
- 4147: LabwareMainTypeDoesNotExist
- 4148: ExperimentDoesNotExist
- 4149: LabwareOrExperimentDoesNotExist
- 4150: CannotAddErrorInfo
- 4151: CurrentLabwareDoesNotMatchPreviousLabware
- 4152: CurrentLabwareAlreadyUsed
- 4153: LabwareIsAlreadyPartOfExperiment
- 4154: ExperimentAlreadyExists
- 4155: DatabaseAlreadyExists
- 4156: DatabaseFileAlreadyAttached
- 4157: CannotFindAValueForAdditionalDataKey
- 4158: IfExistsErrorNotAllowed
- 4159: BadLabwareAccessName
- 4160: CannotFindDeck
- 4161: FileDoesNotExist
- 4162: CannotRemoveExperimentSourceLabware
- 4163: LabwareIsCurrentlyLoadedAndMustBeUnloadedBeforeReloading
- 4164: CannotLookupAdditionalData
- 4165: BarcodeAlreadyUsedAsUniqueBarcode
- 4166: BarcodeNotUnique
- 4167: CannotClearUniqueBarcodeList
- 8193: CannotFindReportDirectory
- 40961: AdditionalDataForeignNotExistsAction
- 40962: AdditionalDataForeignNotExistsInstrument
- 40963: AdditionalDataForeignNotExistsInstrumentConfiguration
- 40964: AdditionalDataForeignNotExistsJob
- 40965: AdditionalDataForeignNotExistsLabware
- 40966: AdditionalDataForeignNotExistsRunAction
- 40967: AdditionalDataForeignNotExistsRun
- 65535: CannotCreateErrorInfo

HxDatabase:

- 0: Unexpected
- 1: CreateTableNotSupported

- 2: CreateProcedureNotSupported
- 3: UnknownEnumValue
- 4: UnsupportedEnumValue
- 5: CannotSynchronizeInternalParameterCollection
- 6: InternalParameterCollectionOutOfSync
- 7: UnknownDBMSIdentifier
- 8: CannotLoadConfigFile
- 9: CannotConvertParameter
- 10: CannotCompareCultureAware
- 11: CannotInferValueFromValue
- 12: BadParameterSupplied
- 13: CannotUseSuppliedConfigFile
- 14: ConfigFileContainsErrors
- 15: ParametersMissing
- 16: SQLScriptCommandTextNotSet
- 17: SQLScriptDoesNotExist
- 18: ExecuteSQLScriptRequiresCommandTypeText
- 19: CannotReadSQLScript
- 20: CannotExecuteSQLScript
- 21: NotImplemented
- 22: CouldNotParseConnectionString
- 65535: CannotCreateErrorInfo

ASWGLOBAL

<https://github.com/theonetrue nerd/VenusPackages/blob/main/ASWGGlobal.pkg>

The ASWGGlobal library doesn't add any functions, but instead declares common HSL constants used in other libraries, SMTs, and methods. These are often used as the return values for functions. The constants it declares are:

- False, No, and Off are all assigned as being `hslFalse`
- True, Yes, and On are all assigned as being `hslTrue`
- **For dialogue buttons, each one is assigned a specific number:**
 - Ok = 1
 - Cancel = 2
 - Abort = 3
 - Retry = 4
 - Ignore = 5
 - Yes = 6
 - No = 7
- Clicking STOP on a timer is assigned as being 3

ALIQOT FOR EASY STEPS V2

https://github.com/theonetrue nerd/VenusPackages/blob/main/Aliquot_for_EasySteps_v2.pkg

This library adds two functions aimed at making aliquoting steps easier. The functions added are:

- `CalcAliquot_v2_1()`
- `CalcChannelPattern()`

CalcAliquot_v2_1(*sequence DispenseSequence, variable VolumePerWell, variable VolumePreAliquote, variable VolumePostAliquote, variable MaxVolumeTip, variable NumberOfDispense, variable NumberOfChannels, variable intNumberOfChannelsInstalled, array arrVolumeToAspirateByChannel, array arrVolumeToMixBeforeAspiration, array arrChannelPatternDispense, variable strChannelPatternAspirate*)

This function calculates the important variables for performing an aliquote pipetting with Easy- or Single Steps

Params DispenseSequence

The sequence for the aliquots to be dispensed into

Params VolumePerWell

The volume of the aliquotes per well

Params VolumePreAliquote

The volume of the pre-aliquote

Params VolumePostAliquote

The volume of the post-aliquote

Params MaxVolumeTip

The maximum volume of the selected CO-RE Tip

Params NumberOfDispense

The number of dispense steps per aspiration

Params NumberOfChannels

The number of the channels to use for the process

Params intNumberOfChannelsInstalled

The number of channels installed in the system

Params arrVolumeToAspirateByChannel

An output of the array of volumes to aspirate by channel

Params arrVolumeToMixBeforeAspiration

An output of the array of volumes to mix before aspiration (typically used during the first aspiration only)

Params arrChannelPatternDispense

The output array of channel patterns for the dispense

Params strChannelPatternAspirate

An output string of the channel pattern to aspirate

Returns

None

Return type

N/A

CalcChannelPattern(*variable i_int_NumberOfPositions, variable i_int_NumberOfInstalledChannels, variable o_str_ChannelPattern*)

This function creates a channel pattern based on the number of desired positions and the number of installed channels.

Params i_int_NumberOfPositions

The number of positions to use, which must be less than or equal to the Number of Installed Channels

Params i_int_NumberOfInstalledChannels

The number of channels installed on the instrument

Params o_str_ChannelPattern

The output channel pattern as a string

Returns

None

Return type

N/A

ALIQOT WITH EASY STEPS

https://github.com/theonetrue nerd/VenusPackages/blob/main/Aliquot_for_EasySteps.pkg

The Aliquot with Easy Steps library adds one function aimed at making aliquoting with Easy- or Single steps easier. This library is outdated, the `Aliquot_for_EasySteps_v2` library should be used instead. The function it adds is:

- `CalcAliquote()`

CalcAliquote(*sequence DispenseSequence, variable VolumePerWell, variable VolumePreAliquote, variable VolumePostAliquote, variable MaxVolumeTip, variable VolumeToAspirate, variable NumberOfDispense, variable FullTrace, variable NumberOfChannels*)

This function calculates the important variables for performing aliquoting with Easy- or Single steps.

Params DispenseSequence

The sequence for the aliquots to be dispensed into

Params VolumePerWell

The volume of the aliquots per well

Params VolumePreAliquote

The volume of the pre-aliquot

Params VolumePostAliquote

The volume of the post-aliquot

Params MaxVolumeTip

The maximum volume of the selected CO-RE Tip

Params VolumeToAspirate

The volume which has to be aspirated

Params NumberOfDispense

The number of dispense steps per aspiration

Params FullTrace

The detail level of the trace (1 is full, 0 is normal)

Params NumberOfChannels

The number of channels installed on the system

Returns

None

Return type

N/A

ARRAY (FROM HSLEXTENSIONS)

<https://github.com/theonetrue nerd/VenusPackages/blob/main/Array.pkg>

The array library from HSLExtensions adds functions to help manipulate 1-D arrays. The following functions are added:

- *Append()*
- *CompareArrays()*
- *Concat()*
- *ContainsDuplicates()*
- *ContainsValue()*
- *ConvertToBooleanArray()*
- *ConvertToFloatArray()*
- *ConvertToIntegerArray()*
- *ConvertToStringArray()*
- *Copy()*
- *FindValue()*
- *InitializeAllValues()*
- *IsBooleanArray()*
- *IsEmpty()*
- *IsFloatArray()*
- *IsIntegerArray()*
- *IsStringArray()*
- *Sort()*

Append(array *io_arrValuesA*, array *i_arrValuesB*)

This function updates the array *io_arrValuesA* to add all the values from *i_arrValuesB* at the end of the array.

Params *io_arrValuesA*

The array to which the values will be added

Params *i_arrValuesB*

The array from which the values will be added

Returns

None

Return type

N/A

CompareArrays(*array i_arrExpectedValues*, *array i_arrActualValues*, *array o_arrMissingValues*, *array o_arrNotExpectedValues*)

This function compares two arrays and outputs arrays of values which are missing from the first array but present in the second, and values which are present in the second array but not in the first.

Params i_arrExpectedValues

The first array, which the second array will be checked against, usually is the array of expected values

Params i_arrActualValues

The second array, which will use the first array as a template when comparing against, usually is your “actual” array

Params o_arrMissingValues

An output array of values which are present in the first array but not the second

Params o_arrNotExpectedValues

An output array of values which are present in the second array but not the first (i.e. unexpected values in your actual data)

Returns

True if both arrays contain the same values (resulting in empty output arrays), false if arrays don’t contain the same values (in which case the output arrays will have data in them)

Return type

Boolean

Concat(*array i_arrValuesA*, *array i_arrValuesB*)

This function appends one array to the other and then returns the concatenated array. The difference between this and the [Append\(\)](#) function is that the Append function updates an existing array, whereas this function doesn’t change the existing arrays and instead returns a new array.

Params i_arrValuesA

The array to which the values will be added

Params i_arrValuesB

The array from which the values will be added

Returns

A new array which is the concatenated version of the input arrays

Return type

Array

ContainsDuplicates(*array i_arrValues*)

This function checks whether the input array has multiple of the same value in it

Params i_arrValues

The array to be checked

Returns

An array with all values which appear more than once in the input array

Return type

Array

ContainsValue(*array i_arrValues*, *variable i_varValue*)

This function determines whether a value exists in an array without returning its index

Params i_arrValues

The array to be searched

Params i_varValue

The value to be searched for

Returns

True if the value is present, false otherwise

Return type

Boolean

ConvertToBooleanArray(array *i_arrValues*, variable *o_blnSuccessfullyConverted*)

This function converts the input array to an array with boolean values. If it is not possible to convert one or more values of the input array, the output will be false and the output array will be empty. Cannot interact with strings, will convert a non-zero int or float into a 1, and will turn a 0 float into a 0.

Params i_arrValues

The array to be converted

Params o_blnSuccessfullyConverted

A boolean which tells you whether the conversion was successful or not

Returns

The boolean version of the input array

Return type

Array

ConvertToFloatArray(array *i_arrValues*, variable *o_blnSuccessfullyConverted*)

This function converts the input array to an array with float values. If it is not possible to convert one or more values of the input array, the output will be false and the output array will be empty. Cannot interact with strings, will convert any int into a float.

Params i_arrValues

The array to be converted

Params o_blnSuccessfullyConverted

A boolean which tells you whether the conversion was successful or not

Returns

The float version of the input array

Return type

Array

ConvertToIntArray(array *i_arrValues*, variable *o_blnSuccessfullyConverted*)

This function converts the input array to one with integer values. If it is not possible to convert one or more values of the input array, the output will be false and the output array will be empty. Cannot interact with strings, will round any floats to the nearest integer.

Params i_arrValues

The array to be converted

Params o_blnSuccessfullyConverted

A boolean which tells you whether the conversion was successful or not

Returns

The integer version of the input array

Return type

Array

ConvertToStringArray(*array i_arrValues*)

This function converts the input array to one with string values.

Params i_arrValues

The array to be converted

Returns

The float version of the input array

Return type

Array

Copy(*array i_arrValues*)

This function will output an exact copy of the input array.

Params i_arrValues

The array to be copied

Returns

A copy of the input array

Return type

Array

FindValue(*array i_arrValues, variable i_varValue*)

This function will lookup the input variable within the input array and return a 1-based array of the indices of all positions that the input variable was found.

Params i_arrValues

The array to be searched

Params i_varValue

The variable to be searched for

Returns

An array of all the locations that the input variable was found

Return type

Array

InitializeAllValues(*array io_arrValues, variable i_varValue*)

This function sets all values within an array to the input variable. Does not work on empty arrays.

Params io_arrValues

The array in which all values will be converted.

Params i_varValue

The variable to which all values will be converted.

Returns

None

Return type

N/A

IsBooleanArray(*array i_arrValues*)

This function checks whether all values in the array are booleans.

Params i_arrValues

The array to be checked

Returns

A boolean of whether the input array is all booleans or not

Return type

Boolean

IsEmpty(array *i_arrValues*)

This function checks whether the input array is empty

Params i_arrValues

The array to be checked

Returns

A boolean of whether the input array is empty or not

Return type

Boolean

IsFloatArray(array *i_arrValues*)

This function checks whether all values in the array are floats.

Params i_arrValues

The array to be checked

Returns

A boolean of whether the input array is all floats or not

Return type

Boolean

IsIntegerArray(array *i_arrValues*)

This function checks whether all values in the array are integers.

Params i_arrValues

The array to be checked

Returns

A boolean of whether the input array is all integers or not

Return type

Boolean

IsStringArray(array *i_arrValues*)

This function checks whether all values in the array are strings.

Params i_arrValues

The array to be checked

Returns

A boolean of whether the input array is all strings or not

Return type

Boolean

Sort(array *i_arrValues*, variable *i_intSortMode*, o_ *bSuccessfulSorted*)

This function outputs a sorted version of the array using the Shakersort sorting algorithm. All values in the array must share the same type for this function to work. Sort mode can either be 1 or 2, 1 is ascending and 2 is descending.

Params i_arrValues

The array containing the values to be sorted

Params i_intSortMode

Whether the array is to be sorted in ascending (1) or descending (2) order

Params o_bSuccessfulSorted

A boolean of whether the sort was successful or not

Returns

A sorted copy of the array

Return type

Array

ARRAYTOOLS

<https://github.com/theonetrue nerd/VenusPackages/blob/main/ArrayTools.pkg>

The ArrayTools library provides the following functions:

- *ArraySeqVLookup()*
- *ArrayVLookup()*
- *ConvertArrayOfNumericIntegersToString()*
- *ConvertArrayOfNumericStringsToInteger()*
- *Lookup()*
- *Sort3ArraysByNumericAscendingOrder()*
- *Update_Value_in_Array()*
- *get_distinct_from_array()*
- *mergeArrays()*
- *removeValueFromArray_basedOnIndex()*

ArraySeqVLookup(*array i_arraySequencesA*, *array i_arrayValuesB*, *variable i_varCondition*, *array o_arraySequencesC*)

Given 2 arrays A (Sequences) and B (Values) with same size Retrieve all the sequences in A where B = input value A(seq1, seq2, seq3, seq4, seq5, seq6, ...) B(1, 2, 3, 1, 2, 3, ...)

Get from A all elements where B = 1 C(seq1, seq4, ...)

Parameters

- **i_arraySequencesA** (*Array of sequences*) – The input array of sequences to filter
- **i_arrayValuesB** (*Array of variables*) – The input array of variables to act as a filter
- **i_varCondition** (*Variable*) – The input value of the variable which is to be selected
- **o_arraySequencesC** (*Array of sequences*) – The output array of the sequences whose index matches the indices in which the the value of the variable is equal to the filter

Returns

None

Return type

N/A

ArrayVLookup(*array i_arrayValuesA, array i_arrayValuesB, variable i_varCondition, array o_arrayValuesC*)

Given 2 arrays of values A and B with same size Retrieve all the elements in B where A = input value

A(100, 20, 15, 45, 42, 1, ...) B(1, 2, 3, 1, 2, 3, ...)

Get from A all elements where B = 1 C(100, 15, ...)

Parameters

- **i_arrayValuesA** (*Array of variables*) – The input array of variables to filter
- **i_arrayValuesB** (*Array of variables*) – The input array of variables to act as a filter
- **i_varCondition** (*Variable*) – The input value of the variable which is to be selected
- **o_arrayValuesC** (*Array of variables*) – The output array of the variables whose index matches the indices in which the value of the variable is equal to the filter

Returns

None

Return type

N/A

ConvertArrayOfNumericIntegersToString(*array i_arr1_int, array o_arr1_str*)

Converts all the numeric integers within an array to strings.

Parameters

- **i_arr1_int** (*Array of variables*) – The input array containing the integers to be converted
- **o_arr1_str** (*Array of variables*) – The output array containing the strings

Returns

None

Return type

N/A

ConvertArrayOfNumericStringsToIntegers(*array i_arr1_str, array o_arr1_int*)

Converts all the numeric strings within an array to integers.

Parameters

- **i_arr1_str** (*Array of variables*) – The input array containing the strings to be converted
- **o_arr1_int** (*Array of variables*) – The output array containing the integers

Returns

None

Return type

N/A

Lookup(*array array, variable item*)

Looks up a value within an array, outputting a 1-based index of the value if found in the array, and a 0 if the value isn't found.

Parameters

- **array** (*Array*) – The input array to be searched
- **item** (*Variable*) – The variable to be searched for

Returns

None

Return type

N/A

Sort3ArraysByNumericAscendingOrder(*array io_array1, array io_array2, array io_array3*)

Sorts 3 arrays by numeric ascending order. *io_array1* must contain only numeric values; this one will be sorted and then the other arrays will update to match the new order of *io_array1*.

Parameters

- **io_array1** (*Array*) – The first of the arrays to be sorted, which must contain only numeric values.
- **io_array2** (*Array*) – The second of the arrays to be sorted, which can contain any values.
- **io_array3** (*Array*) – The third of the arrays to be sorted, which can contain any values.

Returns

None

Return type

N/A

Update_Value_in_Array(*array i_array, variable i_value, variable i_index*)

Overwrites a value in the array at a specified index.

Parameters

- **i_array** (*Array*) – The array in which the value will be changed.
- **i_value** (*Variable*) – The new value to be inserted into the array.
- **i_index** (*Variable*) – The 1-based index of the position in the array to be overwritten.

Returns

None

Return type

N/A

get_distinct_from_array(*array i_arr, array o_arr*)

Gets all the values in an array that only appear once.

Parameters

- **i_arr** (*Array*) – The input array to be searched.
- **o_arr** – The new output array containing the values which only appear once.

Returns

None

Return type

N/A

mergeArrays(*array array1, array array2, array array3*)

Concatenates two arrays and outputs the result into a third array.

Parameters

- **array1** (*Array*) – The first array of interest.
- **array2** (*Array*) – The second array of interest.

- **array3** (*Array*) – The resulting array of values.

Returns

None

Return type

N/A

removeValueFromArray_basedOnIndex(*array i_array_elements, variable i_index_to_remove*)

Removes a value from an array at the specified index.

Parameters

- **i_array_elements** (*Array*) – The array from which an item is to be removed.
- **i_index_to_remove** (*Variable*) – The 1-based index of the array from which the item is to be removed.

Returns

None

Return type

N/A

COLOURSCHEMELIBRARY

<https://github.com/theonetrueerd/VenusPackages/blob/main/ColourSchemeLibrary.pkg>

The ColourSchemeLibrary allows the user to change the colour scheme of the Hamilton standard dialogues. It adds the following function

- *Assign()*

Assign(variable *i_str_ColorSchemeCode*)

This function updates the colour scheme of the Standard Dialogues that will be used in the method where the library is executed. Requires HSLExtensions(File), HSLExtensions(Directory), and ASW Global.

Params i_str_ColorSchemeCode

The input color scheme code which is the int KitProviderCode. Hamilton = 0.

Returns

None

Return type

N/A

CONVERT FILE TO ASCII

<https://github.com/theonetrueerd/VenusPackages/blob/main/ConvertFileToASCII.pkg>

The Convert File to ASCII library adds the ability to convert Hamilton files to ASCII files. The function it adds is:

- *ConvertFileToASCII()*

ConvertFileToASCII(*variable i_Str_SourceFile*)

This function converts files from binary to ASCII. The file must be a Hamilton file.

Params i_Str_SourceFile

The path to the file which is being converted.

Returns

None

Return type

N/A

DEBUGASSIST

<https://github.com/theonetrue nerd/VenusPackages/blob/main/DebugAssist.pkg>

The debug assist library adds one function aimed at helping the user identify, understand, and fix runtime errors. The function it adds is:

- *FullErrorAnalysis()*

This function can use the [ASWStandardDialog library](#) which makes the dialog nicer, but is not a requirement.

FullErrorAnalysis(*variable i_PathForErrorDataAnalysis*)

This function is designed to be called in the “OnAbort” submethod. It will look at the error that caused the abort to be triggered, convert the trace file error code into the more standard form, identify what that error corresponds with and hopefully suggest some initial things to check. In order for the dialogue to pop up, this function requires ASWStandardDialogues to be initialised.

Params i_PathForErrorDataAnalysis

A file path for where the error data will be exported. Currently not very important as the main bonus of the library is the dialogue that pops up, although the intention is to add more detail to this exported file.

Returns

None

Return type

N/A

EDITFILEATTRIBUTES

<https://github.com/theonetrue nerd/VenusPackages/blob/main/EditFileAttributes.pkg>

The Edit File Attributes library adds the following functions:

- *Make_Hidden()*
- *Make_ReadOnly()*
- *Remove_Hidden()*
- *Remove_ReadOnly()*

Make_Hidden(*variable Filename*)

Applies the hidden attribute to the file found at the path specified.

Params filename

The path to the file that is going to be modified.

Returns

None

Return type

N/A

Make_ReadOnly(*variable Filename*)

Applies the read only attribute to the file found at the path specified.

Params filename

The path to the file that is going to be modified.

Returns

None

Return type

N/A

Remove_Hidden(*variable Filename*)

Removes the hidden attribute from the file found at the path specified.

Params filename

The path to the file that is going to be modified.

Returns

None

Return type

N/A

Remove_ReadOnly(*variable Filename*)

Removes the read only attribute from the file found at the path specified.

Params filename

The path to the file that is going to be modified.

Returns

None

Return type

N/A

ERROR REPORT LIBRARY

<https://github.com/theonetrue nerd/VenusPackages/blob/main/ErrorReportLibrary.pkg>

The error report library generates concise reports on both general errors and pipetting errors that occur during a run, extracting the information from the trace file. It adds the following functions:

- `CopyActiveTraceFile()`
- `GenerateGeneralErrorsReport()`
- `GeneratePipettingErrorsReport()`

CopyActiveTraceFile(*variable o_strTraceCopyFilePath*)

Finds the trace file that is currently active (i.e. associated with the current run) and makes a copy of it, which it saves in the location specified.

Params o_strTraceCopyFilePath

The location that the copy of the trace file will be created in

Returns

None

Return type

N/A

GenerateGeneralErrorsReport(*variable i_strTraceFilePath, variable o_strPDF_FilePath*)

Reads the active trace file, locates errors found within it and generates an error report on errors found in non-pipetting steps

Params i_strTraceFilePath

The path to the current trace file

Params o_strPDF_FilePath

The file path for the exported PDF of the error report

Returns

None

Return type

N/A

GeneratePipettingErrorsReport(*variable i_strTraceFilePath, variable i_strTemplateFilePath, variable i_intWriteRowStart, o_strPDF_FilePath*)

Reads the active trace file, locates errors in the pipetting steps within it and generates an error report for the pipetting steps only

Params i_strTraceFilePath

The path to the current trace file

Params i_strTemplateFilePath

The path to the excel file used as a template for generating the report; should be installed when the library is unpacked

Params i_intWriteRowStart

Which row of the excel file to begin the report on; usually two, can be more than two

Params o_strPDF_FilePath

The file path for the exported PDF of the error report

Returns

None

Return type

N/A

ERROR SIMULATOR

<https://github.com/theonetrue nerd/VenusPackages/blob/main/ErrorSimulator.pkg>

This library adds the ability to simulate a variety of errors, to assist with debugging or custom error handling. The functions it adds are:

- `AA_Abstract()`
- `STEP1_PrepareRegistryAndCfgFile()`
- `STEP2a_SimulateError_Channels()`
- `STEP2b_SimulateError_COREGripper()`
- `STEP2c_SimulateError_iSWAP()`
- `STEP2d_SimulateError_96Head()`
- `STEP2e_SimulateError_384Head()`
- `STEP2f_SimulateError_BarcodeReading()`
- `STEP2g_SimulateError_Autoload()`
- `STEP2h_SimulateError_CRWashstation()`
- `STEP3_Restore_BackupCfgFile()`
- `STEP4_Optional_SwitchChecksum_ON()`

AA_Abstract()

This function does nothing. In the submethod library itself it simply lists the changelog.

STEP1_PrepareRegistryAndCfgFile(device ML_STAR)

This function prepares the system files and registry to perform the error simulation. You can put this step at the beginning of your method, and you can disable it after the 1st run, once all registry and configuration settings are done. It switches off the checksum in the registry, asking for confirmation to modify the registry during run. It converts the ConfigML_STAR_Simulator.cfg file to ASCII. It creates a copy of ML_STAR_Simulator.cfg named ML_STAR_Simulator.cfg.bak for future restoration if needed. It replaces the string reload “0” with reload “1”.

Params ML_STAR

The instrument being used. Should be ML_STAR.

Returns

None

Return type

N/A

STEP2a_SimulateError_Channels(*variable ChannelNumber, variable WhenSimulateError, variable ErrorToSimulate, device ML_STAR*)

This function simulates an error in the channels. It can simulate an error of your choice (between errors 1 and 12) in a desired channel at any step of the pipetting process. The error codes are: - 1..... Liquid Level not found (error 06/70) - 2..... Not enough liquid (error 06/71) - 3..... Liquid level not found with Dual LLD (error 06/73) - 4..... Clot error (error 04/81) - 5..... Liquid not correctly aspirated (error 06/80) - 6..... Tip already present (error 07/00) - 7..... No Tip (error 08/00) - 8..... Wrong tip type detected (error 08/78) - 9..... Syntax / parameter out of range (error 01/00) - 10..... Hardware (error 02/00) - 11..... Not Executed (error 03/00) - 12..... Not Completed (error 10/00)

Params ChannelNumber

A string of the channel(s) in which the error is to be simulated, e.g. "1,2,4,7"

Params WhenSimulateError

The step of the process in which the error is desired to be simulated. 1 = Aspiration, 2 = Dispense, 3 = Tip pickup, 4 = Tip eject

Params ErrorToSimulate

Integer corresponding to the desired simulated error, key listed in function description.

Params ML_STAR

The instrument being used. Should be ML_STAR.

Returns

None

Return type

N/A

STEP2b_SimulateError_COREGripper(*variable ChannelNumber, variable WhenSimulateError, variable ErrorToSimulate, device ML_STAR*)

This function simulates an error in the CO-RE gripper. It can simulate an error of your choice (between errors 1 and 4) in a desired channel at any step of the process. The error codes are: - 1..... Can't get the CORE grippers (error 08/75) - 2..... Step lost in Z drive, crash against something (error 02/62) - 3..... Hardware (error 02/00) - 4..... Read Barcode Error (05/00)

Params ChannelNumber

A string of the channel(s) in which the error is to be simulated, e.g. "1,2,4,7"

Params WhenSimulateError

The step of the process in which the error is desired to be simulated. 1 = Pickup CO-RE, 2 = Get Plate, 3 = Put Plate, 4 = Eject CO-RE gripper

Params ErrorToSimulate

Integer corresponding to the desired simulated error, key listed in function description.

Params ML_STAR

The instrument being used. Should be ML_STAR.

Returns

None

Return type

N/A

STEP2c_SimulateError_iSWAP(*variable WhenSimulateError, variable ErrorToSimulate, device ML_STAR*)

This function simulates an error in the iSWAP. It can simulate an error of your choice (between errors 1 and 4) at any step of the process. The error codes are: - 1..... Expected object not found (error 21/94) - 2..... Step lost in Z drive, crash against something (error 02/62) - 3..... Hardware (error 02/00) - 4..... Object lost (error 23/96)

Params WhenSimulateError

The step of the process in which the error is desired to be simulated. 1 = Get Plate, 2 = Put Plate

Params ErrorToSimulate

Integer corresponding to the desired simulated error, key listed in function description.

Params ML_STAR

The instrument being used. Should be ML_STAR.

Returns

None

Return type

N/A

STEP2d_SimulateError_96Head(*variable WhenSimulateError, variable ErrorToSimulate, device ML_STAR*)

This function simulates an error in the 96-head. It can simulate an error of your choice (between errors 1 and 12 - not all available) at any step of the process. The error codes are: - 1..... Liquid Level not found (error 06/70) - 2..... Not enough liquid (error 06/71) - 6..... Tip already present (error 07/00) - 7..... No Tip (error 08/00) - 8..... Wrong tip type detected (error 08/78) - 9..... Syntax / parameter out of range (error 01/00) - 10..... Hardware (error 02/00) - 11..... Not Executed (error 03/00) - 12..... Not Completed (error 10/00)

Params WhenSimulateError

The step of the process in which the error is desired to be simulated. 1 = Aspiration, 2 = Dispense, 3 = Tip Pickup, 4 = Tip eject, 5 = Wash

Params ErrorToSimulate

Integer corresponding to the desired simulated error, key listed in function description.

Params ML_STAR

The instrument being used. Should be ML_STAR.

Returns

None

Return type

N/A

STEP2e_SimulateError_384head(*variable WhenSimulateError, variable ErrorToSimulate, device ML_STAR*)

This function simulates an error in the 384-head. It can simulate an error of your choice (between errors 1 and 12 - not all available) at any step of the process. The error codes are: - 1..... Liquid Level not found (error 06/70) - 2..... Not enough liquid (error 06/71) - 6..... Tip already present (error 07/00) - 7..... No Tip (error 08/00) - 8..... Wrong tip type detected (error 08/78) - 9..... Syntax / parameter out of range (error 01/00) - 10..... Hardware (error 02/00) - 11..... Not Executed (error 03/00) - 12..... Not Completed (error 10/00)

Params WhenSimulateError

The step of the process in which desired error is to be simulated. 1 = Aspiration, 2 = Dispense, 3 = Tip Pickup, 4 = Tip eject, 5 = Wash.

Params ErrorToSimulate

Integer corresponding to the desired simulated error, key listed in function definition.

Params ML_STAR

The instrument being used. Should be ML_STAR

Returns

None

Return type

N/A

STEP2f_SimulateError_BarcodeReading(*variable WhenSimulateError, device ML_STAR*)

This function simulates an error during the barcode reading step. It can simulate the error as though it was either the CO-RE grips or the iSWAP involved.

Params WhenSimulateError

Poorly labelled variable, likely copied and pasted. This is not when the error is simulated; it is instead whether the barcode is being read with CO-RE grips (1) or the iSWAP (2)

Parameters

ML_STAR (*Device*) – The instrument being used. Should be ML_STAR

Returns

None

Return type

N/A

STEP2g_SimulateError_Autoload(*variable notReadPosition_Str, variable notPresentPosition_Str, device ML_STAR*)

This submethod simulates the “Barcode not read” and “Labware not Present Errors”.you can enter positions with errors separated by a comma

Params notReadPosition_Str

A string of comma separated values with the positions in which you want to simulate barcode errors

Params notPresentPosition_Str

A string of comma separated values with the positions in which you want to simulate labware not present errors.

Params ML_STAR

The instrument being used. Should be ML_STAR.

Returns

None

Return type

N/A

STEP2h_SimulateError_CRWashstation(*variable errorToSimulate, device ML_STAR*)

This function simulates an error when starting the washing with the CR washstation

Params errorToSimulate

The error being simulated. 1 = Not washing liquid error, 2 = hardware error.

Params ML_STAR

The instrument being used. Should be ML_STAR.

Returns

None

Return type

N/A

STEP3_Restore_BackupCfgFile(*device ML_STAR*)

This function restores the backup configuration file created in function STEP1_PrepRegAndCfgFile

Params ML_STAR

The instrument being used. Should be ML_STAR.

Returns

None

Return type

N/A

STEP4_Optional_SwitchChecksum_ON()

This function is optional. Use it at the end of your method if you want to Switch the File Checksum ON in the registry

Returns

None

Return type

N/A

FRAMEWORK (FROM HSLEXTENSIONS)

<https://github.com/theonetrue nerd/VenusPackages/blob/main/Framework.pkg>

The framework library from HSLExtensions provides two visible functions, but is also required for a few of the other HSLExtension libraries. It uses the functionality of the ASW Standard TraceLevel library. The functions it adds are:

- *GetVersion()*
- *SetTraceLevel()*

GetVersion()

This function gets the framework version

Returns

The framework version

Return type

String

SetTraceLevel(*variable i_intTraceLevel*)

This function sets the trace level of the framework.

Params i_intTraceLevel

The trace level, for the ASWStandard::TraceLevel library. 0 = TRACE_LEVEL_NONE, 1 = TRACE_LEVEL_RELEASE, 2 = TRACE_LEVEL_DEBUG.

Returns

None

Return type

N/A

HSLDECKVISUALIZE

<https://github.com/theonetrueerd/VenusPackages/blob/main/HSLDeckVisualize.pkg>

The HSLDeckVisualize library adds functions which will update the layout shown in the run control to account for manually performed actions or cause specific effects to occur. The functions it adds are:

- *UpdateUsedPositions()*
- *UpdateUsedLabware()*
- *UpdateLoadedLabware()*

UpdateUsedPositions(*device dev, sequence positions, variable action, variable description*)

This function updates the list of used positions and updates the view.

Params dev

The instrument for which the associated deck view will be updated. ML_STAR should be the only option visible in the dropdown.

Params positions

The sequence of labware and positions IDs to be updated.

Params action

The action to be performed on the specified positions. The same action will be performed for all positions. Action can be an integer from 0-6, with the values each representing a different event. 0 = Position is selected, 1 = Processing, 2 = Reserved, 3 = Error, 4 = Processed, 5 = Reset action state to none, 6 = Position selected and flashing.

Params description

The description displayed in the first section of the status bar in the deck view.

Returns

None

Return type

N/A

UpdateUsedLabware(*device dev, array labware, array action, variable description*)

This function updates the list of used labware and updates the view accordingly. The labware and action arrays must be the same size.

Params dev

The instrument for which the associated deck view will be updated. ML_STAR should be the only option visible in the dropdown.

Params labware

The array of labware IDs for updating

Params action

The array of actions which will be applied to the associated labware at the same index in the labware array. Action can be an integer from 0-6, with the values each representing a different event. 0 = Position is selected, 1 = Processing, 2 = Reserved, 3 = Error, 4 = Processed, 5 = Reset action state to none, 6 = Position selected and flashing.

Params description

The description displayed in the first section of the status bar in the deck view.

Returns

None

Return type

N/A

UpdateLoadedLabware(*device dev, array labware, array state, variable description*)

This function updates the list of loaded labware and updates the view accordingly. The labware and action arrays must be the same size.

Params dev

The instrument for which the associated deck view will be updated. ML_STAR should be the only option visible in the dropdown.

Params labware

The array of labware IDs for updating

Params state

The array of load state which will be applied to the associated labware at the same index in the labware array. State can be an integer from 0-5, with the values each representing a different state. 0 = unload the labware and make not visible, 1 = load the labware, 2 = preparing to unload, 3 = preparing to load, 4 = labware will flash, 5 = labware will flash

Params description

The description displayed in the first section of the status bar in the deck view.

Returns

None

Return type

N/A

HSLFILLIB

<https://github.com/theonetrue nerd/VenusPackages/blob/main/HSLFilLib.pkg>

This library allows interaction with and manipulation of files present on the host computer. It adds the following functions:

- *FileOf()*
- *FilFindFile()*
- *FilFindNextFile()*
- *FilFormatBarcodeFile()*
- *FilGetBinPath()*
- *FilGetCommState()*
- *FilGetCommTimeouts()*
- *FilGetConfigPath()*
- *FilGetLabwarePath()*
- *FilGetLibraryPath()*
- *FilGetLogFilesPath()*
- *FilGetMethodsPath()*
- *FilGetSystemPath()*
- *FilIsNull()*
- *FilReadString()*
- *FilRemoveFields()*
- *FilSearchPath()*
- *FilSetCommState()*
- *FilSetCommTimeouts()*
- *FilUpdateRecord()*
- *FilWriteString()*

FileOf(*variable filObj*)

This function checks whether the current position in the specified file is the final line

Params filObj

The opened file to be checked

Returns

Boolean as to whether the position is the end of the file or not

Return type

Boolean

FilFindFile(*variable fileName*)

This function starts searching the specified path for a file. This function is obsolete, and [FilSearchPath\(\)](#) should be used instead.

Params fileName

The directory or path and file name to be searched for. Can contain wildcard characters such as * or ?

Returns

If successful, the path name of the first file found that matches the input

Return type

Variable

FilFindNextFile()

Continues the search from [FilFindFile\(\)](#) to the next file.

Returns

If successful, the path name of the next file found that matches the input from the most recent FilFindFile command

Return type

Variable

FilFormatBarcodeFile(*variable dataSource, variable dataTarget*)

This function takes the barcode ASCII text file written during LoadCarrier and converts it into a strongly formatted barcode file. This strongly formatted file can be an ASCII text file, a Microsoft Excel file, or a Microsoft Access file. Will contain the following columns:

- ID (record ID, integer)
- Specifier (string, P = position, C = carrier)
- Position (position, integer)
- Barcode (barcode, string)
- Timestamp (timestamp, YYYY-MM-DD hh:mm:ss, string)

Params dataSource

The name of the barcode ASCII file generated during the load carrier step

Params dataTarget

The name of the barcode ASCII text file, Microsoft Excel file or Microsoft Access file generated by the function. The name must include a table name for a Microsoft Excel file or a Microsoft Access Database.

Returns

Boolean showing if the function was successful or not

Return type

Boolean

FilGetBinPath()

This function retrieves the vector binary path

Returns

The vectory binary path (usually C:Program Files (x86)HamiltonBin)

Return type

Variable

FilGetCommState(*file port*)

This function retrieves the configuration information for the specified communication resource. The entries of the structure that retrieves the configuration information must be accessible in the global scope.

Params port

The communication resource opened during the file-Open operation

Type

Port

Returns

Boolean showing if the function was successful or not

Return type

Boolean

GetCommTimeouts(*file port*)

This function retrieves the time-out parameters for all read and write operations for the specified communication resource. The entries of the structure that contains the configuration information must be accessible in the global scope.

Params port

The communication resource opened during the file-Open operation

Returns

Boolean showing if the function was successful or not

Return type

Boolean

FilGetConfigPath()

This function retrieves the vector configuration path

Returns

The vector configuration path (usually C:Program Files (x86)HamiltonConfig)

Return type

Variable

FilGetLabwarePath()

This function retrieves the vector labware path

Returns

The vector labware path (usually C:Program Files (x86)HamiltonLabware)

Return type

Variable

FilGetLibraryPath()

This function retrieves the vector library path

Returns

The vector library path (usually C:Program Files (x86)HamiltonLibrary)

Return type

Variable

FilGetLogFilePath()

This function retrieves the vector log files path

Returns

The vector log files path (usually C:Program Files (x86)HamiltonLogFiles)

Return type

Variable

FilGetMethodsPath()

This function retrieves the vector methods path

Returns

The vector methods path (usually C:Program Files (x86)HamiltonMethods)

Return type

Variable

FilGetSystemPath()

This function retrieves the vector system path

Returns

The vector system path (usually C:Program Files (x86)HamiltonSystem)

Return type

Variable

FilIsNull(*variable value*)

This function returns a non-zero if the variable is a null value (SQL style Null).

Params value

The variable being checked

Returns

A boolean determining if the variable is SQL style Null or not

Return type

Boolean

FilReadString(*file fileObj*)

This function reads the next record from the input file as string-valued data. Row data, but no schema data, is saved to the string. After you call FilReadString, the next unread record becomes the current record, or Eof is set to hslTrue if there are no more records.

Params fileObj

The file being looked at

Returns

The contents of the line being looked at as string-valued data, or the specific run-time error

Return type

String

FilRemoveFields(*file fileObj*)

This function removes all fields from a record definition

Params fileObj

The file containing the record which is having its fields removed

Returns

None

Return type

N/A

FilSearchPath(*variable fileName*)

This function searches for the specified file, and outputs the path and filename of the file if found, or an empty string if not found. Will search the current directory, the methods directory, the library directory, and any directories in the PATH environment variable.

Params fileName

The file name to be searched for

Returns

The path name of the first file found, or an empty string if no files were found

Return type

Variable

FilSetCommState(*file port, variable cfgFile*)

This function configures a communication resource according to the specifications in a structure that contains the configuration information. The structure that contains the configuration information must be structured as shown below. Each entry in the structure is optional and overwrites the default value in parentheses.

Params port

The communication resource opened during the file-Open operation

Params cfgFile

The name of the file containing the configuration information. If this parameter is empty, the entries in the structure that contains the configuration information must be accessible in the global scope.

Returns

Boolean showing whether the function succeeded or not

Return type

Boolean

FilSetCommTimeouts(*file port, variable cfgFile*)

This function sets the time-out parameters for all read and write operations on a specified communication resource. The structure that contains the time-out information is as shown below. Each entry in the structure is optional and overwrites the default value in parentheses.

Params port

The communication resource opened during the file-Open operation

Params cfgFile

The name of the file that contains the time-out information. If this parameter is empty, the entries in the structure that contains the time-out information must be accessible in the global scope.

Returns

Boolean showing whether the function succeeded or not

Return type

Boolean

FilUpdateRecord(*file fileObj*)

Updates the current record of the file object with the values of the variable objects specified in the record definition. The current record remains current after you call the FilUpdateRecord function. The provider must support UPDATE.

Params fileObj

The file object being updated

Returns

Boolean showing whether the function succeeded or not

Return type

Boolean

FilWriteString(*file fileObj, variable stringObj*)

Writes a string to the end of the file data source. After you call the FilWriteString function, the new record becomes the current record.

Params fileObj

The file which is being written in

Params stringObj

The string to be written

Returns

Boolean showing whether the function succeeded or not

Return type

Boolean

HSLFILELIBEX

<https://github.com/theonetrueerd/VenusPackages/blob/main/HSLFilLibEx.pkg>

HSLFilLibEx allows you to manipulate files. It adds the following functions:

- *FilCopyFileEx()*
- *FilDeleteFileEx()*
- *FilFormatReportFileEx()*

FilCopyFileEx(*variable SourceFilePathName, variable DestinationFilePathName*)

This function creates a copy of the specified file, which will be placed in the specified new location

Params SourceFilePathName

The path to the file which you wish to copy, including the file name and extension

Params DestinationFilePathName

The path to which you would like to copy the file, including the file name and extension

Returns

None

Return type

N/A

FilDeleteFileEx(*variable FilePathName*)

This function deletes the file at the specified location

Params FilePathName

The path to the file which you wish to delete, including file name and extension

Returns

None

Return type

N/A

FilFormatReportFileEx(*variable SourceFilePathName, variable DestinationFilePathName*)

This function creates a formatted/filtered copy of the specified report file. It searches the file for the string “Element Name” and copies any subsequent lines to the new file.

Params SourceFilePathName

The path to the file which you wish to format, including file name and extension

Params DestinationFilePathName

The path to the file you would like to create the formatted report in, including file name and extension.

Returns

None

Return type

N/A

HSLSTATISTICS

<https://github.com/theonetrueerd/VenusPackages/blob/main/HSLStatistics.pkg>

The HSLStatistics library is designed to easily allow the user to perform basic statistical functions easily. It adds the following functions:

- *Stat_Average()*
- *Stat_StdDeviation()*
- *Stat_CorrelationCoefficient()*
- *Stat_RSQ()*
- *Stat_Slope()*
- *Stat_Intercept()*
- *Stat_Min()*
- *Stat_Max()*

Stat_Average(array DataArray)

This function returns the average of a data set input as an array

Params DataArray

The input array of data from which the average is calculated

Returns

The calculated average of the array

Return type

Variable

Stat_StdDeviation(array DataArray)

This function returns the standard deviation of a data set input as an array

Params DataArray

The input array of data from which the standard deviation is calculated

Returns

The calculated standard deviation of the array

Return type

Variable

Stat_CorrelationCoefficient(array XArray, array YArray)

This function returns the correlation coefficient (r-value) of a paired data set

Params XArray

The array of X data values in the paired set

Params YArray

The array of Y data values in the paired set

Returns

The r-value of the paired data set

Return type

Variable

Stat_RSQ(*array XArray, array YArray*)

This function returns the pearson coefficient (r^2) of a paired data set

Params XArray

The array of X data values in the paired set

Params YArray

The array of Y data values in the paired set

Returns

The pearson coefficient of the paired data set

Return type

Variable

Stat_Slope(*array XArray, array YArray*)

This function returns the slope of the best fit line of a paired data set

Params XArray

The array of X data values in the paired set

Params YArray

The array of Y data values in the paired set

Returns

The slope of the best fit line

Return type

Variable

Stat_Intercept(*array XArray, array YArray*)

This function returns the intercept of the best fit line of a paired data set

Params XArray

The array of X data values in the paired set

Params YArray

The array of Y data values in the paired set

Returns

The intercept of the line of best fit

Return type

Variable

Stat_Min(*array DataArray*)

This function returns the lowest value of a dataset

Params DataArray

The array for the minimum to be searched in

Returns

The lowest value of the dataset

Return type
Variable

HSLSTRLIB

<https://github.com/theonetrueenerd/VenusPackages/blob/main/HSLStrLib.pkg>

HSL String Library provides the following functions:

- *StrAsciiToStr()*
- *StrConcat2()*
- *StrConcat4()*
- *StrConcat8()*
- *StrConcat12()*
- *StrEvaluateExpr()*
- *StrFillLeft()*
- *StrFillRight()*
- *StrFind()*
- *StrFStr()*
- *StrFStrEx()*
- *StrFVal()*
- *StrGetLength()*
- *StrGetType()*
- *StrHexIStr()*
- *StrIsDigit()*
- *StrIStr()*
- *StrIVal()*
- *StrLeft()*
- *StrMakeLower()*
- *StrMakeLowerCopy()*
- *StrMakeUpper()*
- *StrMakeUpperCopy()*
- *StrMid()*
- *StrReplace()*

- *StrReverseFind()*
- *StrRight()*
- *StrSpanExcluding()*
- *StrStrToAscii()*
- *StrTrimLeft()*
- *StrTrimRight()*

StrAsciiToStr(*variable asciiCode*)

Converts the given ASCII Code (an integer) to a character (string).

Parameters

asciiCode (*Integer*) – The ASCII code to convert

Returns

The ASCII code as a string

Return type

String

StrConcat2(*variable Var1, variable Var2*)

Combines two strings into a new string

Parameters

- **Var1** (*String*) – The first string to be combined
- **Var2** (*String*) – The second string to be combined

Returns

The combined string of Var1 + Var2

Return type

String

StrConcat4(*variable Var1, variable Var2, variable Var3, variable Var4*)

Combines four strings into a new string

Parameters

- **Var1** (*String*) – The first string to be combined
- **Var2** (*String*) – The second string to be combined
- **Var3** (*String*) – The third string to be combined
- **Var4** (*String*) – The fourth string to be combined

Returns

The combined string of Var1 + Var2 + Var3 + Var4

Return type

String

StrConcat8(*variable Var1, variable Var2, variable Var3, variable Var4, variable Var5, variable Var6, variable Var7, variable Var8*)

Combines eight strings into a new string

Parameters

- **Var1** (*Variable*) – The first string to be combined

- **Var2** (*Variable*) – The second string to be combined
- **Var3** (*Variable*) – The third string to be combined
- **Var4** (*Variable*) – The fourth string to be combined
- **Var5** (*Variable*) – The fifth string to be combined
- **Var6** (*Variable*) – The sixth string to be combined
- **Var7** (*Variable*) – The seventh string to be combined
- **Var8** (*Variable*) – The eighth string to be combined

Returns

The combined string of Var1 + Var2 + Var3 + Var4 + Var5 + Var6 + Var7 + Var8

Return type

String

StrConcat12(*variable Var1, variable Var2, variable Var3, variable Var4, variable Var5, variable Var6, variable Var7, variable Var8, variable Var9, variable var10, variable var11, variable var12*)

Combines twelve strings into a new string

Parameters

- **Var1** (*Variable*) – The first string to be combined
- **Var2** (*Variable*) – The second string to be combined
- **Var3** (*Variable*) – The third string to be combined
- **Var4** (*Variable*) – The fourth string to be combined
- **Var5** (*Variable*) – The fifth string to be combined
- **Var6** (*Variable*) – The sixth string to be combined
- **Var7** (*Variable*) – The seventh string to be combined
- **Var8** (*Variable*) – The eighth string to be combined
- **Var9** (*Variable*) – The ninth string to be combined
- **Var10** (*Variable*) – The tenth string to be combined
- **Var11** (*Variable*) – The eleventh string to be combined
- **Var12** (*Variable*) – The twelfth string to be combined

Returns

The combined string of Var1 + Var2 + Var3 + Var4 + Var5 + Var6 + Var7 + Var8 + Var9 + Var10 + Var11 + Var12

Return type

Variable

StrEvaluateExpr(*variable expression*)

This function evaluates an expression within a string. All variables involved must have global scope.

Params expression

The expression to evaluate as a string

Returns

The value of the expression if the function succeeds, otherwise a runtime error

Return type

Variable

StrFillLeft(*variable str, variable character, variable width*)

This function fills leading characters to the string

Params str

The string to be modified

Params character

The user-defined character to be filled

Params width

The width to be filled

Returns

The modified string

Return type

Variable

StrFillRight(*variable str, variable character, variable width*)

This function fills trailing characters to the string

Params str

The string to be modified

Params character

The user-defined character to be filled

Params width

The width to be filled

Returns

The modified string

Return type

Variable

StrFind(*variable str, variable subStr*)

This function searches the string for the first match of the sub-string.

Params str

The string to be searched

Params subStr

The substring to be searched for

Returns

The zero-based index of the first character in this string object that matches the requested sub-string or characters. -1 if the sub-string is not found.

Return type

Integer

StrFStr(*variable number*)

This function converts the floating point number input into the corresponding character string.

Params number

The float to be converted into a string

Returns

The string form of the float

Return type

Variable

StrFStrEx(*variable number*, *variable languageSpecific*, *variable precision*)

This function converts the floating point number input into the corresponding character string

Params number

The float to be converted into a string.

Params languageSpecific

Boolean which specifies whether the decimal symbol in the Regional Settings should be used to write the string representation of the floating point number.

Params precision

The total number of significant digits to be used for the floating-point display

Returns

The string representation of the floating point number

Return type

Variable

StrFVal(*variable str*)

Converts the sequence of digits, contained in the character string str, into the corresponding floating point number. Conversion aborts at the first character in str, which is not a digit or not one of the characters +, -, e, E.

Params str

The string to be converted

Returns

The float representation of the input string. Null if the string cannot be converted. DBL_MAX if the conversion results in an overflow. DBL_MIN if the conversion results in an underflow.

Return type

Float or variable

StrGetLength(*variable str*)

Returns the number of characters in a string object (without '0').

Params str

The string being read

Returns

The length of the string

Return type

Integer

StrGetType(*variable var*)

This function retrieves the type of the value of a variable

Params var

A reference to a variable (int, float or string)

Returns

One of the following string-valued constants that indicates the type of the value of a variable. i = hslInteger, f = hslFloat, s = hslString, null = no type

Return type

Variable

StrHexIStr(*variable number*)

Converts the input integer into the corresponding hexadecimal character string

Params number

The integer to be converted

Type

Integer

Returns

The hexadecimal string representation of the integer

Return type

String

StrIsDigit(*variable character*)

The StrIsDigit function determines if the specified input string (which should be a character) is a digit or not.

Params character

The input character to be a tested, as a string

Returns

Boolean showing whether the character is a digit (1) or not (0)

Return type

Boolean

StrIsStr(*variable number*)

The StrIsStr function converts the input integer into the corresponding character string

Params number

The integer to be converted

Returns

The string representation of the integer number

Return type

Variable

StrIVal(*variable str*)

The StrIVal function converts the input sequence of digits into the corresponding integer. The input string is treated as a decimal, unless it begins with an 0x in which case it is interpreted as hexadecimal. Conversion aborts at the first character in the input which is neither a digit nor one of the characters “+” or “-”.

Params str

The input sequence of digits to be converted

Returns

The numeric value of the sequence of digits contained in the character string, as an integer. Null if the character string cannot be converted into a number. LONG_MAX = 2147483647 if the conversion results in an overflow. LONG_MIN = -2147483647 - 1 if the conversion results in an underflow.

Return type

Variable

StrLeft(*variable str, variable count*)

The StrLeft function extracts the first (leftmost) characters of a string and returns a copy of the extracted substring. The number of characters extracted is equal to the input variable “count”. If “count” is longer than the string, the entire string is returned.

Params str

The input string from which the substring is to be extracted

Params count

The number of characters to be extracted

Returns

A string containing a copy of the specified range of characters. Can be an empty string.

Return type

Variable

StrMakeLower(*variable str*)

The StrMakeLower function converts the original string to its lowercase form.

Params str

The string to be converted

Returns

The original string converted to lowercase

Return type

Variable

StrMakeLowerCopy(*variable str*)

The StrMakeLowerCopy function returns a copy of the original string converted to lowercase.

Params str

The string to be copied

Returns

A copy of the original string converted to lowercase

Return type

Variable

StrMakeUpper(*variable str*)

The StrMakeUpper function converts the original string to its uppercase form.

Params str

The string to be converted

Returns

The original string converted to uppercase

Return type

Variable

StrMakeUpperCopy(*variable str*)

The StrMakeUpperCopy function returns a copy of the original string converted to uppercase.

Params str

The string to be copied

Returns

A copy of the original string converted to uppercase

Return type

Variable

StrMid(*variable str, variable first, variable count*)

The StrMid function extracts a substring of length “count” characters from the input variable “str”, starting at position “first” which is 0-based. The function returns a copy of the extracted substring.

Params str

The string from which the substring is to be extracted

Params first

The first character to be extracted (0-based)

Params count

The number of characters to be extracted

Returns

A string containing a copy of the specified range of characters, can be empty

Return type

Variable

StrReplace(*variable str, variable oldSubStr, variable newSubStr*)

The StrReplace function searches a string for a specified substring and replaces it with another specified substring.

Params str

The string to be edited

Params oldSubStr

The substring to be replaced by newSubStr

Params newSubStr

The substring to replace oldSubStr

Returns

The number of replaced instances of oldSubStr. Zero if the string is unchanged.

Return type

Variable

StrReverseFind(*variable str, variable subStr*)

The StrReverseFind function searches a string object for the last match of a sub-string

Params str

The string to be searched

Params subStr

The substring to be searched for

Returns

The zero-based index of the last character in this string that matches the requested substring or characters. -1 if the substring is not found.

Return type

Variable

StrRight(*variable str, variable count*)

The StrRight function extracts the last (rightmost) characters of a string and returns a copy of the extracted substring. The number of characters extracted is equal to the input variable “count”. If “count” is longer than the string, the entire string is returned.

Params str

The input string from which the substring is to be extracted

Params count

The number of characters to be extracted

Returns

A string containing a copy of the specified range of characters. Can be an empty string.

Return type

Variable

StrSpanExcluding(*variable str, variable subStr*)

The StrSpanExcluding function can be used to search the string for the first occurrence of any character in the specified set subStr. StrSpanExcluding extracts and returns all characters preceding the first occurrence of a character from subStr (in other words, the character from subStr and all characters following it in the string, are not returned). If no character from subStr is found in the string, then StrSpanExcluding returns the entire string.

Params str

The string to be searched

Params subStr

A string containing the set of characters to be searched for

Returns

A sub-string containing characters in the string that are not in subStr, beginning with the first character in the string and ending with the first character found in the string that is also in subStr (that is, starting with the first character in the string and up to but excluding the first character in the string that is found subStr). It returns the entire string if no character in subStr is found in the string.

Return type

Variable

StrStrToAscii(*variable character*)

The StrStrToAscii function converts the given character (as a string) into an ASCII code (as an integer)

Params character

The character to convert, inputted as a string

Returns

The ASCII code for the given character as an integer, -1 if the function fails

Return type

Variable

StrTrimLeft(*variable str, variable character*)

The StrTrimLeft function trims leading whitespace characters from the string (removes newline, space, tab, and user-defined characters)

Params str

The string to trim

Params character

A string containing user-defined characters to be trimmed (may be empty, in which case only newline, space and tabs will be trimmed)

Returns

None

Return type

N/A

StrTrimRight(*variable str, variable character*)

The StrTrimRight function trims lagging whitespace characters from the string (removes newline, space, tab, and user-defined characters)

Params str

The string to trim

Params character

A string containing user-defined characters to be trimmed (may be empty, in which case only newline, space and tabs will be trimmed)

Returns

None

Return type

N/A

HSLZIPLIB

<https://github.com/theonetrue nerd/VenusPackages/blob/main/HSLZipLib.pkg>

HSLZipLib is a library that allows you to create zip folders, as well as unzip and extract folders from them. It adds the following four functions:

- *UnpackArchive()*
- *AddItemToArchive()*
- *InitializeArchive()*
- *PackArchive()*

UnpackArchive(*variable i_strArchiveName, variable i_strPathOutput*)

The unpack archive function unzips a zip folder and extracts the contents to the specified location.

Params i_strArchiveName

The full name and path to the zip folder you wish to unzip, including the .zip file extension

Params i_strPathOutput

The full name and path of the folder you would like to create when extracting the zip file, with no extension

Returns

Boolean confirming whether unzipping the archive was successful or not

Return type

Boolean

AddItemToArchive(*variable i_strFileOrDirectoryName, variable i_strDirectoryPathInArchive*)

This function specifies a file or directory which will be added to the zip file which will be generated with the *PackArchive()* function.

Params i_strFileOrDirectoryName

The path to the file or directory that you wish to include in the zip folder. File extension required for files, no extension required for folders.

Params i_strDirectoryPathInArchive

The path to the desired location of the file within the zip folder.

Returns

Boolean confirming whether adding the file or directory to the archive was successful

Return type

Boolean

InitializeArchive(*variable i_strArchiveName*)

This function initializes whatever archive is specified. It is required to call this function for each archive you are interacting with, in advance of calling any of the otehr functions.

Params i_strArchiveName

The name of the desired zip file; can be a file that already exists or one that doesn't exist yet and will be created

Returns

None

Return type

N/A

PackArchive()

This function zips/packs whichever archive was most recently initialized.

Returns

Boolean confirming whether zipping the archive was successful or not

Return type

Boolean

HSL_SEQDAILYTOOLS

https://github.com/theonetrueerd/VenusPackages/blob/main/HSL_SeqDailyTools.pkg

The HSL_SeqDailyTools library adds four functions aimed at making some parts of sequence handling slightly easier. The four functions it adds are:

- `CopyPlatePattern96ToTipRack()`
- `CopyPlatePatternToPlate()`
- `GetNumberOfPositionsLeft()`
- `SeqEasyEdit()`

CopyPlatePattern96ToTipRack(*sequence plateSource, sequence tipRack*)

This function copies the sequence pattern on a 96 well plate and applies it to a tip rack. It is used to ensure that the 96 head only picks up tips in the wells that the plate is going to be interacting with.

Params plateSource

The sequence pattern on the plate to be copied

Params tipRack

The sequence of the tip rack which will have the template applied to it

Returns

None

Return type

N/A

CopyPlatePatternToPlate(*sequence plateSource, sequence plateTarget*)

This function copies the pattern of wells from one plate to another plate. It is useful when a plate is being moved around on deck, or when reagents are being pipetted from one plate into the same locations on another plate.

Params plateSource

The sequence pattern from the plate to be copied

Params plateTarget

The sequence of the second plate which is having the template applied to it

Returns

None

Return type

N/A

GetNumberOfPositionsLeft(*sequence seq, variable numberOfPositionsLeft*)

This function calculates the number of positions left in a sequence

Params seq

The sequence which is being checked

Params numberOfPositionsLeft

The output value of how many positions are left in the sequence

Returns

None

Return type

N/A

SeqEasyEdit (*sequence seq, variable timeout, variable dialog_title, variable dialog_msg, device ML_STAR*)

This function is a simplified version of SeqEdit, with only the parameters that are regularly used. It opens a dialog box from which the user can edit whichever sequence is specified in the input.

Params seq

The sequence which is being edited

Params timeout

How long the dialog window will stay open by itself before it closes and continues the method without editing the sequence

Params dialog_title

The title of the dialog window which pops up

Params dialog_message

The text within the dialog window which pops up

Params ML_STAR

The device used. Should be ML_STAR from the dropdown

Returns

None

Return type

N/A

IF_AND_IF_OR

https://github.com/theonetrue nerd/VenusPackages/blob/main/If_And_If_Or.pkg

The `if_and_if_or` library adds two functions:

- `if_and()`
- `if_or()`

`if_and`(*variable iVar1, variable iVar2, variable iVar3, variable iVar4*)

This function checks variable `iVar1` against `iVar2` and `iVar3` against `iVar4`; if both statements are true then the function returns a 1, otherwise it returns a 0

Params iVar1

Variable to be checked against `iVar2`

Params iVar2

Variable to be checked against `iVar1`

Params iVar3

Variable to be checked against `iVar4`

Params iVar4

Variable to be checked against `iVar3`

Returns

1 if both true, 0 otherwise

Return type

Variable

`if_or`(*variable iVar1, variable iVar2, variable iVar3, variable iVar4*)

This function checks variable `iVar1` against `iVar2` and `iVar3` against `iVar4`; if at least one of the statements is true then the function will return a 1, otherwise it will return a 0

Params iVar1

Variable to be checked against `iVar2`

Params iVar2

Variable to be checked against `iVar1`

Params iVar3

Variable to be checked against `iVar4`

Params iVar4

Variable to be checked against `iVar3`

Returns

1 if at least one is true, 0 if false

Return type
Variable

LABWARE PROPERTIES

https://github.com/theonetrueerd/VenusPackages/blob/main/Labware_Properties.pkg

The labware properties library adds a variety of functions which assist with obtaining the physical data of the labware, as well as things like its ID. The functions it adds are:

- `GetCarrierIDandSiteID_FromLabID()`
- `Get_ContainerBaseOffset()`
- `Get_ContainerBaseThickness()`
- `Get_Height()`
- `Get_NameAndFileName()`
- `Get_NumberOfColumns()`
- `Get_NumberOfRows()`
- `Get_RackBaseToCoverBase()`
- `Get_StackHeight()`
- `Get_XYZ_deckPosition()`
- `Get_XY_dimensions()`

GetCarrierIDandSiteID_FromLabID(*device io_instrument, variable i_labware_ID, variable o_carrier_ID, variable o_site_ID*)

This function takes outputs the carrier ID and site ID which are associated with the given labware ID. It works on both the Nimbus and the STAR.

Params io_instrument

The instrument being used. Will either be ML_STAR or Nimbus.

Params i_labware_ID

The labware ID from which the other IDs are being pulled.

Params o_carrier_ID

The carrier ID associated with the given labware ID.

Params o_site_ID

The site ID associated with the given labware ID.

Returns

A boolean determining whether the labware ID exists (1) or not (0)

Return type

Boolean

Get_ContainerBaseOffset(*device io_instrument, sequence i_sequenceLabware, variable i_sequencePosition, variable o_ContainerBaseOffset*)

This function outputs the distance from the rack base to the container base for the labware at the given sequence position.

Params io_instrument

The instrument being used. Will either be ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware being checked.

Params i_sequencePosition

The position ID or index of the sequence being checked. Can be either int or str. An input of 0 will auto-select the first position.

Params o_ContainerBaseOffset

The distance from the rack base to the container base at the specified position

Returns

None

Return type

N/A

Get_ContainerBaseThickness(*device io_instrument, sequence i_sequenceLabware, variable o_containerBaseThickness*)

This function outputs the base thickness of the container at the first position of a given sequence.

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware being checked.

Params o_containerBaseThickness

The thickness of the base of the container being checked.

Returns

None

Return type

N/A

Get_Height(*device io_instrument, sequence i_sequenceLabware, variable o_labwareHeight*)

This function outputs the height of the labware at the first position of a given sequence. This value is only the labware height, not the absolute Z position.

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware being checked.

Params o_labwareHeight

The height of the labware being checked.

Returns

None

Return type

N/A

Get_NameAndFileName(*device io_instrument, sequence i_sequenceLabware, variable o_viewName, variable o_fileName*)

This function outputs the labware view name and the file name associated with it (with path)

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware of interest.

Params o_viewName

The view name of the labware being checked.

Params o_fileName

The file name (with path) of the labware being checked.

Returns

None

Return type

N/A

Get_NumberOfColumns(*device io_instrument, sequence i_sequenceLabware, variable o_labwareColumns*)

This function outputs the number of columns defined in the labware from the first position of a given sequence.

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware of interest.

Params o_labwareColumns

The number of columns defined in the labware being checked.

Returns

None

Return type

N/A

Get_NumberOfRows(*device io_instrument, sequence i_sequenceLabware, variable o_labwareColumns*)

This function outputs the number of rows defined in the labware at the first position of a given sequence. The variable and function description both say columns; this is incorrect.

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware being checked.

Params o_labwareColumns

The number of rows defined in the labware being checked.

Returns

None

Return type

N/A

Get_RackBaseToCoverBase(*device io_instrument, sequence i_sequenceLabware, variable o_RackBaseToCoverBase_Height*)

This function outputs the height from the base of the rack to the base of the cover/lid of the labware at the first position of a given sequence.

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware being checked.

Params o_RackBaseToCoverBase_Height

The distance from the rack base to the cover base.

Returns

None

Return type

N/A

Get_StackHeight(*device io_instrument, sequence i_sequenceLabware, variable o_labwareStackHeight*)

This function outputs the stack height of the specified labware at the first position of a given sequence.

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware being checked.

Params o_labwareStackHeight

The stack height of the labware being checked, or the covered stack height if the labware is lidded.

Returns

None

Return type

N/A

Get_XYZ_deckPosition(*device io_instrument, sequence i_sequenceLabware, variable o_labware_deckPosition_X, variable o_labware_deckPosition_Y, variable o_labware_deckPosition_Z*)

This function returns the X, Y and Z coordinates of the upper left well of the specified labware at the first position of a given sequence. The description of this function says it only does the X and Y coordinates, this is incorrect.

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware being checked.

Params o_labware_deckPosition_X

The X coordinate of the labware on the deck.

Params o_labware_deckPosition_Y

The Y coordinate of the labware on the deck.

Params o_labware_deckPosition_Z

The Z coordinate of the labware on the deck.

Returns

None

Return type

N/A

Get_XY_dimensions(*device io_instrument, sequence i_sequenceLabware, variable o_X_width, variable o_Y_depth*)

This function outputs the X (width) and Y (depth) dimensions of the specified labware at the first position of a given sequence.

Params io_instrument

The instrument being used. Will be either ML_STAR or Nimbus.

Params i_sequenceLabware

The sequence associated with the labware being checked.

Params o_X_width

The width of the labware being checked.

Params o_Y_depth

The depth of the labware being checked.

Returns

None

Return type

N/A

LOOKUP

<https://github.com/theonetrue nerd/VenusPackages/blob/main/Lookup.pkg>

The lookup library provides one function which allows you to locate items in an array and get their index.

- *Lookup()*

Lookup(*array array, variable item*)

Input an array and a value to look up. If the value occurs in the array, it will return the 1-based index. If it doesn't occur in the array, it will return a 0.

Params array

The array to search

Params item

The variable to search for

Returns

1-based index if found, 0 if not found

Return type

Variable

PIPETTING (FROM HSLEXTENSIONS)

<https://github.com/theonetrueerd/VenusPackages/blob/main/Pipetting.pkg>

The pipetting library from HSL Extensions adds functions related to pipetting, primarily to do with sequences and channel patterns. The functions it adds are:

- *CreateMixedChannelPattern()*
- *CreateOrderedChannelPattern()*
- *ParseChannelPattern()*
- *Plate384PositionNumberToPositionText()*
- *Plate384PositionTextToPositionNumber()*
- *Plate96PositionNumberToPositionText()*
- *Plate96PositionTextToPositionNumber()*

CreateMixedChannelPattern(array *i_arrUseChannel*, variable *i_intTotalNumberOfChannels*)

This function generates a “mixed” channel pattern. The channel pattern generated will be based on the input array, and then will be filled with “0”s up to the total number of channels (normally 8). If the array size is greater than the total number of channels inputted, the channel pattern will be limited to whatever the total number of channels inputted is. E.g. if the array is [1,1,1,1,0,0,1,0,0,1] and the input channel number is 8, the channel pattern generated will be “11110010”.

Params *i_arrUseChannel*

A boolean array of which channels should be used.

Params *i_intTotalNumberOfChannels*

The total number of channels being used.

Returns

The channel pattern generated by the function

Return type

String

CreateOrderedChannelPattern(variable *i_intNumberOfUsedChannels*, variable *i_intTotalNumberOfChannels*)

This function creates an ordered channel pattern, generating a string with as many “1”s in it as the inputted number of used channels, and adding “0”s until the total length of the channel pattern is equal to the total number of channels inputted. If the total number of used channels is greater than the total number of channels, the channel pattern will be truncated to the total number of channels. E.g. if the number of used channels is 3 and the total number of channels is 8, the channel pattern generated will be “11100000”. If these numbers were reversed, the channel pattern generated would simply be “111”. :params *i_intNumberOfUsedChannels*: The number of “active” channels wishing to be used in the channel pattern :params *i_intTotalNumberOfChannels*: The max number of channels available (usually 8 on a STAR) :type *i_intNumberOfUsedChannels*: Variable :type *i_intTotalNumberOfChannels*: Variable :return: The channel pattern generated by the function :rtype: String

ParseChannelPattern(*variable i_strChannelPattern*)

This function parses a channel pattern given as a string and creates a boolean array of the channels being used. If the input has the wrong type, or contains not allowed characters, the result is an empty array.

Params i_strChannelPattern

The channel pattern to be parsed

Returns

The array of booleans with the channels being used in the channel pattern

Return type

Array

Plate384PositionNumberToPositionText(*variable i_intPositionNumber, variable o_strPositionText*)

This function converts position number to position text for a 384 well plate. E.g. if given a 1 it will convert it to A1, if given a 16 it will convert it to P1, etc.

Params i_intPositionNumber

The position number to be converted. Any int in the range of 1-384.

Params o_strPositionText

The position text which will be the output. Will be a string of any of the well identifiers on a 384 well plate (i.e. A1-P24)

Returns

A boolean of whether the conversion was successful

Return type

Boolean

Plate384PositionTextToPositionNumber(*variable i_strPositionText, variable o_intPositionNumber*)

This function converts position text to position number for a 384 well plate. E.g. if given A1 it will convert it to a 1, if given P1 it will convert it to a 16, etc.

Params i_strPositionText

The position text to be converted. Will be a string of any of the well identifiers on a 384 well plate (i.e. A1-P24)

Params o_intPositionNumber

The position number which will be the output. Any int in the range of 1-384.

Returns

A boolean of whether the conversion was successful

Return type

Boolean

Plate96PositionNumberToPositionText(*variable i_intPositionNumber, variable o_strPositionText*)

This function converts position number to position text for a 96 well plate. E.g. if given a 1 it will convert it to A1, if given an 8 it will convert it to H1, etc.

Params i_intPositionNumber

The position number to be converted. Any int in the range of 1-96.

Params o_strPositionText

The position text which will be the output. Will be a string of any of the well identifiers on a 96 well plate (i.e. A1-H12)

Returns

A boolean of whether the conversion was successful

Return type

Boolean

Plate96PositionNumberToPositionText(*variable i_strPositionText, variable o_intPositionNumber*)

This function converts position text to position number for a 96 well plate. E.g. if given A1 it will convert it to a 1, if given H1 it will convert it to a 8, etc.

Params i_strPositionText

The position text to be converted. Will be a string of any of the well identifiers on a 96 well plate (i.e. A1-H12)

Params o_intPositionNumber

The position number which will be the output. Any int in the range of 1-96.

Returns

A boolean of whether the conversion was successful

Return type

Boolean

STAR_CHANNEL_TOOLS

https://github.com/theonetrueerd/VenusPackages/blob/main/STAR_Channel_Tools.pkg

The STAR_Channel_Tools submethod library adds a variety of functions which are related to the use of the 8 channels. The following functions are added:

- *CHAN_ACCESS_Sort1Sequence()*
- *CHAN_ACCESS_Sort1Sequence1Array()*
- *CHAN_ACCESS_Sort1Sequence2Arrays()*
- *CHAN_ACCESS_Sort2Sequences()*
- *CHAN_ACCESS_Sort2Sequences1Array()*
- *CHAN_ACCESS_Sort2Sequences2Arrays()*
- *LIQUID_LEVEL_GetLiquidLevelHeight()*
- *LIQUID_LEVEL_MeasureLiquidMulti()*
- *LIQUID_LEVEL_MeasureLiquidSingle()*
- *LIQUID_LEVEL_ReturnVolumesFromLiquidLevel()*
- *MOVE_ChannelsToSequencePosition()*
- *MOVE_ChannelsToSequencePosition_5mL()*
- *MOVE_CheckPlateWithTwoChannels()*
- *MOVE_InitDispenseDrive()*
- *MOVE_InitDispenseDrive_5mL()*
- *PLATE_STACK_CountPlateStacks()*
- *QUERY_GetChannelPosition()*
- *QUERY_GetChannelPosition_5mL()*
- *QUERY_GetTipPresentState()*
- *QUERY_GetTipPresentState_5mL()*
- *QUERY_GetTipVolume()*
- *QUERY_GetTipVolume_5mL()*
- *SPLIT_WELLS_AddContainersToWell()*
- *SPLIT_WELLS_RemoveContainers()*
- *TRAVEL_LANES_MoveChannelsToTravelLanes()*

- `TRAVEL_LANES_MoveChannelsToTravelLanes_5mL()`
- `TRAVEL_LANES_MoveChannelsToYPosition()`
- `TRAVEL_LANES_MoveChannelsToYPosition_5mL()`
- `TRAVEL_LANES_MoveChannelsWithTravelLanes()`
- `TRAVEL_LANES_MoveChannelsWithTravelLanes_5mL()`
- `TRAVEL_LANES_SingleSource_ChannelDisplacement()`
- `TRAVEL_LANES_SingleSource_ChannelDisplacement_5mL()`

CHAN_ACCESS_Sort1Sequence(*device ML_STAR, sequence io_Sequence_to_Sort, variable i_Channel_Type, boolean i_Sort_by_Labware, boolean i_Sort_by_XY, boolean i_Sort_for_Channel_Raster, variable i_Max_Channel, sequence o_Sorted_Sequence, variable o_Channel_Pattern*)

This submethod takes an input sequence and sorts it based on the input parameters of labware, position, and raster. Once sorted, the submethod will choose a position that the current channel can access up to the maximum. If the current channel cannot access the position, it will skip it and move to the next available position. If the current cannot access any more positions, that channel will be skipped. Make sure the channel use setting is set to “All Sequence Positions” in the pipetting step, otherwise the sequence and channel pattern will not line up.

Params ML_STAR

The ML_STAR itself, which will be the only option in the dropdown.

Params io_Sequence_to_Sort

The input sequence to be sorted.

Params i_Channel_Type

The channel type associated with the pipetting step (1mL, 5mL, labware handler). 0 = 1mL, 1 = 5mL, 2 = Labware handler.

Params i_Sort_by_Labware

A boolean determining whether the sequence is to be sorted by labware in ascending order

Params i_Sort_by_XY

A boolean determining whether the sequence is to be sorted by position (X ascending, Y descending)

Params i_Sort_for_Channel_Raster

A boolean determining whether the next position will be at least the raster distance unless no other positions are available

Params i_Max_Channel

The maximum channel that you want to be used from 1 to 16. 0 turns this option off and the maximum number of channels will be used.

Params o_Sorted_Sequence

The outputted sorted sequence for the pipetting step

Params o_Channel_Pattern

The outputted channel pattern for the pipetting step

Returns

The number of sequence positions remaining in the sequence

Return type

Variable

CHAN_ACCESS_Sort1Sequence1Array(*device ML_STAR, sequence io_Sequence_to_Sort, array io_Array_of_Variables, variable i_Channel_Type, boolean i_Sort_by_Labware, boolean i_Sort_by_XY, boolean i_Sort_for_Channel_Raster, variable i_Max_Channel, sequence o_Sorted_Sequence, array o_Sorted_Array, variable o_Channel_Pattern*)

This submethod takes in input sequence and sorts it by the conditions given below. After sorting, the submethod will choose a position that the current channel can access up to the maximum. If the current channel cannot access the position, it will skip it and move to the next available position. If the current channel cannot access any more positions, that channel will be skipped. The array will be sorted with the sequence. The array and the sequence must be the same size. Make sure the channel use setting is set to “All sequence positions” otherwise the sequence and channel pattern will not line up.

Params ML_STAR

The ML_STAR itself, which will be the only option in the dropdown.

Params io_Sequence_to_Sort

The input sequence to be sorted.

Params io_Array_of_Variables

The array to be sorted with the sequence. Used positions will be removed from the array.

Params i_Channel_Type

The channel type associated with the pipetting step (1mL, 5mL, labware handler). 0 = 1mL, 1 = 5mL, 2 = Labware handler.

Params i_Sort_by_Labware

A boolean determining whether the sequence is to be sorted by labware in ascending order

Params i_Sort_by_XY

A boolean determining whether the sequence is to be sorted by position (X ascending, Y descending)

Params i_Sort_for_Channel_Raster

A boolean determining whether the next position will be at least the raster distance unless no other positions are available

Params i_Max_Channel

The maximum channel that you want to be used from 1 to 16. 0 turns this option off and the maximum number of channels will be used.

Params o_Sorted_Sequence

The outputted sorted sequence for the pipetting step

Params o_Sorted_Array

The outputted sorted array which matches the sequence

Params o_Channel_Pattern

The outputted channel pattern for the pipetting step

Returns

The number of sequence positions remaining in the sequence

Return type

Variable

CHAN_ACCESS_Sort1Sequence2Arrays(*device ML_STAR, sequence io_Sequence_to_Sort, array io_Array_of_Variables, array io_Array_of_Variables2, variable i_Channel_Type, boolean i_Sort_by_Labware, boolean i_Sort_by_XY, boolean i_Sort_for_Channel_Raster, variable i_Max_Channel, sequence o_Sorted_Sequence, array o_Sorted_Array, array o_Sorted_Array2, variable o_Channel_Pattern*)

This submethod takes in input sequence and sorts it by the conditions given below. After sorting, the submethod will choose a position that the current channel can access up to the maximum. If the current channel cannot access the position, it will skip it and move to the next available position. If the current channel cannot access any more positions, that channel will be skipped. The arrays will be sorted with the sequence. The arrays and the sequence must be the same size. Make sure the channel use setting is set to “All sequence positions” otherwise the sequence and channel pattern will not line up.

Params ML_STAR

The ML_STAR itself, which will be the only option in the dropdown.

Params io_Sequence_to_Sort

The input sequence to be sorted.

Params io_Array_of_Variables

The first array to be sorted with the sequence. Used positions will be removed from the array.

Params io_Array_of_Variables2

The second array to be sorted with the sequence. Used positions will be removed from the array.

Params i_Channel_Type

The channel type associated with the pipetting step (1mL, 5mL, labware handler). 0 = 1mL, 1 = 5mL, 2 = Labware handler.

Params i_Sort_by_Labware

A boolean determining whether the sequence is to be sorted by labware in ascending order

Params i_Sort_by_XY

A boolean determining whether the sequence is to be sorted by position (X ascending, Y descending)

Params i_Sort_for_Channel_Raster

A boolean determining whether the next position will be at least the raster distance unless no other positions are available

Params i_Max_Channel

The maximum channel that you want to be used from 1 to 16. 0 turns this option off and the maximum number of channels will be used.

Params o_Sorted_Sequence

The outputted sorted sequence for the pipetting step

Params o_Sorted_Array

The outputted sorted first array which matches the sequence

Params o_Sorted_Array2

The outputted sorted second array which matches the sequence

Params o_Channel_Pattern

The outputted channel pattern for the pipetting step

Returns

The number of sequence positions remaining in the sequence

Return type

Variable

CHAN_ACCESS_Sort2Sequences(*device ML_STAR, sequence io_Sequence_to_Sort, sequence io_Sequence_to_Sort2, variable i_Channel_Type, boolean i_Sort_by_Labware, boolean i_Sort_by_XY, boolean i_Sort_for_Channel_Raster, variable i_Max_Channel, sequence o_Sorted_Sequence, sequence o_Sorted_Sequence2, variable o_Channel_Pattern*)

This sub method takes the in input sequences and sorts them by the conditions given below. After sorting, the sub will choose a position that the current channel can access in both sequence positions up to the maximum. If the current channel cannot access the positions, it will skip it and move to the next available position. If the current channel cannot access any more positions, that channel will be skipped. Make sure the channel use setting is set to “All sequence positions” otherwise the sequence and channel pattern will not line up. The first sequence is the driving sequence and the second sequence will be adjusted by the first sort.

Params ML_STAR

The ML_STAR itself, which will be the only option in the dropdown.

Params io_Sequence_to_Sort

The first input sequence to be sorted.

Params io_Sequence_to_Sort2

The second input sequence to be sorted.

Params i_Channel_Type

The channel type associated with the pipetting step (1mL, 5mL, labware handler). 0 = 1mL, 1 = 5mL, 2 = Labware handler.

Params i_Sort_by_Labware

A boolean determining whether the sequence is to be sorted by labware in ascending order

Params i_Sort_by_XY

A boolean determining whether the sequence is to be sorted by position (X ascending, Y descending)

Params i_Sort_for_Channel_Raster

A boolean determining whether the next position will be at least the raster distance unless no other positions are available

Params i_Max_Channel

The maximum channel that you want to be used from 1 to 16. 0 turns this option off and the maximum number of channels will be used.

Params o_Sorted_Sequence

The outputted second sorted sequence for the pipetting step

Params o_Sorted_Sequence2

The outputted second sorted sequence for the pipetting step

Params o_Channel_Pattern

The outputted channel pattern for the pipetting step

Returns

The number of sequence positions remaining in the sequence

Return type

Variable

CHAN_ACCESS_Sort2Sequences1Array(device ML_STAR, sequence io_Sequence_to_Sort, sequence io_Sequence_to_Sort2, array io_Array_of_Variables, variable i_Channel_Type, boolean i_Sort_by_Labware, boolean i_Sort_by_XY, boolean i_Sort_for_Channel_Raster, variable i_Max_Channel, sequence o_Sorted_Sequence, sequence o_Sorted_Sequence2, array o_Sorted_Array, variable o_Channel_Pattern)

This sub method takes the in input sequences and sorts them by the conditions given below. After sorting, the sub will choose a position that the current channel can access in both sequence positions up to the maximum. If the current channel cannot access the positions, it will skip it and move to the next available position. If the current channel cannot access any more positions, that channel will be skipped. Make sure the channel use setting is set

to “All sequence positions” otherwise the sequence and channel pattern will not line up. The first sequence is the driving sequence and the second sequence will be adjusted by the first sort.

Params ML_STAR

The ML_STAR itself, which will be the only option in the dropdown.

Params io_Sequence_to_Sort

The first input sequence to be sorted.

Params io_Sequence_to_Sort2

The second input sequence to be sorted.

Params i_Channel_Type

The channel type associated with the pipetting step (1mL, 5mL, labware handler). 0 = 1mL, 1 = 5mL, 2 = Labware handler.

Params i_Sort_by_Labware

A boolean determining whether the sequence is to be sorted by labware in ascending order

Params i_Sort_by_XY

A boolean determining whether the sequence is to be sorted by position (X ascending, Y descending)

Params i_Sort_for_Channel_Raster

A boolean determining whether the next position will be at least the raster distance unless no other positions are available

Params i_Max_Channel

The maximum channel that you want to be used from 1 to 16. 0 turns this option off and the maximum number of channels will be used.

Params o_Sorted_Sequence

The outputted second sorted sequence for the pipetting step

Params o_Sorted_Sequence2

The outputted second sorted sequence for the pipetting step

Params o_Sorted_Array

The sorted array. The array will be the size of the maximum channel.

Params o_Channel_Pattern

The outputted channel pattern for the pipetting step

Returns

The number of sequence positions remaining in the sequence

Return type

Variable

CHAN_ACCESS_Sort2Sequences2Arrays(*device ML_STAR, sequence io_Sequence_to_Sort, sequence io_Sequence_to_Sort2, array io_Array_of_Variables, array io_Array_of_Variables2, variable i_Channel_Type, boolean i_Sort_by_Labware, boolean i_Sort_by_XY, boolean i_Sort_for_Channel_Raster, variable i_Max_Channel, sequence o_Sorted_Sequence, sequence o_Sorted_Sequence2, array o_Sorted_Array, array o_Sorted_Array2, variable o_Channel_Pattern*)

This sub method takes the in input sequences and sorts them by the conditions given below. After sorting, the sub will choose a position that the current channel can access in both sequence positions up to the maximum. If the current channel cannot access the positions, it wil skip it and move to the next available position. If the current channel cannot access any more positions, that channel will be skipped. Make sure the channel use setting is set

to “All sequence positions” otherwise the sequence and channel pattern will not line up. The first sequence is the driving sequence and the second sequence will be adjusted by the first sort.

Params ML_STAR

The ML_STAR itself, which will be the only option in the dropdown.

Params io_Sequence_to_Sort

The first input sequence to be sorted. This positions removed will be removed from this sequence.

Params io_Sequence_to_Sort2

The second input sequence to be sorted. The positions used will be removed from this sequence.

Params io_Array_of_Variables

The input array of variables which will be sorted with the first sequence. The positions used will be removed from the array.

Params io_Array_of_Variables2

The input array of variables which will be sorted with the second sequence. The positions used will be removed from the array.

Params i_Channel_Type

The channel type associated with the pipetting step (1mL, 5mL, labware handler). 0 = 1mL, 1 = 5mL, 2 = Labware handler.

Params i_Sort_by_Labware

A boolean determining whether the sequence is to be sorted by labware in ascending order

Params i_Sort_by_XY

A boolean determining whether the sequence is to be sorted by position (X ascending, Y descending)

Params i_Sort_for_Channel_Raster

A boolean determining whether the next position will be at least the raster distance unless no other positions are available

Params i_Max_Channel

The maximum channel that you want to be used from 1 to 16. 0 turns this option off and the maximum number of channels will be used.

Params o_Sorted_Sequence

The outputted second sorted sequence for the pipetting step

Params o_Sorted_Sequence2

The outputted second sorted sequence for the pipetting step

Params o_Sorted_Array

The sorted array. The array will be the size of the maximum channel.

Params o_Sorted_Array2

The second sorted array. The array will be the size of the maximum channel.

Params o_Channel_Pattern

The outputted channel pattern for the pipetting step

Returns

The number of sequence positions remaining in the sequence

Return type

Variable

LIQUID_LEVEL_GetLiquidLevelHeight(*device ML_STAR, variable i_str_LiquidLevelReturn, sequence i_seq_Labware, variable i_int_Channel, variable o_flt_LiquidHeight*)

This function will return the liquid level height relative to the container bottem.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown.

Params i_str_LiquidLevelReturn

The return value of the liquid level detect from the pipetting step.

Params i_seq_Labware

The input sequence from which the height is to be determined

Params i_int_Channel

The channel which will be used to determine the liquid level height

Params o_flt_LiquidHeight

The detected liquid level height

Returns

None

Return type

N/A

LIQUID_LEVEL_MeasureLiquidMulti(*device ML_STAR, array i_arrseq_FullReservoirSequences, sequence i_seq_TipsToUse, sequence i_seq_TipWaste, variable i_str_TipCounter, variable i_int_LLD_Sensitivity, variable i_bool_ConvertTouL, array o_arr_VolumesMeasured*)

This function will pick up the desired tips and will measure the liquid level at the center most well of the desired reservoirs and will return the volumes in uL. Ensure the sequence used for the reservoir contains the FULL number of positions of the reservoir, otherwise volume estimation will be incomplete! The tip types supported are: 50uL Filter, 50uL, 300uL Filter, 300uL, 1000uL Filter, and 1000uL.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown.

Params i_arrseq_FullReservoirSequences

The array of full reservoir sequences to be checked.

Params i_seq_TipsToUse

The sequence of the tips to be used to measure the liquid level.

Params i_seq_TipWaste

The sequence of the waste used to eject the tips.

Params i_str_TipCounter

The tip counter to be used for the tips on pickup. Place "" if no tip counter will be used.

Params i_int_LLD_Sensitivity

Integer representing the liquid level sensitivity to be used. 1 is very high, 4 is low, 5 is from labware definition.

Params i_bool_ConvertTouL

Boolean determining whether it converts to uL (hs!True) or stays as mL (hs!False).

Params o_arr_VolumesMeasured

The array of volumes that were measured.

Returns

Whether the measurement has been successful or not

Return type

Boolean

LIQUID_LEVEL_MeasureLiquidSingle(*device ML_STAR, sequence i_seq_FullReservoirSequence, sequence i_seq_TipsToUse, sequence i_seq_TipWaste, variable i_str_TipCounter, variable i_bool_IncrementTipSequence, variable i_int_LLD_Sensitivity, variable i_bool_ConvertTouL, variable o_flt_VolumeMeasured*)

This function will pick up the desired tip and will measure the liquid level at the center most well of the desired reservoir and will return the volume in uL. Ensure the sequence used for the reservoir contains the FULL number of positions of the reservoir, otherwise volume estimation will be incomplete! The tip types supported are: 50uL Filter, 50uL, 300uL Filter, 300uL, 1000uL Filter, and 1000uL.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown.

Params i_seq_FullReservoirSequences

The full reservoir sequence to be checked.

Params i_seq_TipsToUse

The sequence of the tips to be used to measure the liquid level.

Params i_seq_TipWaste

The sequence of the waste used to eject the tips.

Params i_str_TipCounter

The tip counter to be used for the tips on pickup. Place "" if no tip counter will be used.

Params i_bool_IncrementTipSequence

Boolean determining whether the tip sequence should be incremented after pickup or not.

Params i_int_LLD_Sensitivity

Integer representing the liquid level sensitivity to be used. 1 is very high, 4 is low, 5 is from labware definition.

Params i_bool_ConvertTouL

Boolean determining whether it converts to uL (hslTrue) or stays as mL (hslFalse).

Params o_flt_VolumeMeasured

A float of the volume that was measured.

Returns

Whether the measurement has been successful or not

Return type

Boolean

LIQUID_LEVEL_ReturnVolumesFromLiquidLevel(*device ML_STAR, variable i_str_PipettingReturn, variable i_str_LiquidLevelReturn, variable i_bool_ConvertTouL, array o_arr_VolumesMeasured*)

This function will return the volumes that were measured from a previous aspiration step.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_str_PipettingReturn

The return value of the pipetting step to measure liquid level

Params i_str_LiquidLevelReturn

The return value of the liquid level detect

Params i_bool_ConvertTouL

Boolean determining whether it converts to uL (hslTrue) or stays as mL (hslFalse)

Params o_arr_VolumesMeasured

An array of the volumes which were measured for each channel

Returns

None

Return type

N/A

MOVE_ChannelsToSequencePosition(*device ML_STAR, variable i_str_ChPattern, sequence i_seq_Positions, variable i_flt_ZHeight*)

This function moves the 1mL channels to set positions. This function will only move the channels that are activated by the channel pattern. The positions in the sequence will skip over the positions where the channel is turned off.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_str_ChPattern

The channel pattern of the pipettes to move

Params i_seq_Positions

The X Y positions to move the channels to

Params i_flt_ZHeight

The Z positions to end the channels in

Returns

None

Return type

N/A

MOVE_ChannelsToSequencePosition_5mL(*device ML_STAR, variable i_str_ChPattern, sequence i_seq_Positions, variable i_flt_ZHeight*)

This function moves the 5mL channels to set positions. This function will only move the channels that are activated by the channel pattern. The positions in the sequence will skip over the positions where the channel is turned off.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_str_ChPattern

The channel pattern of the pipettes to move

Params i_seq_Positions

The X Y positions to move the channels to

Params i_flt_ZHeight

The Z positions to end the channels in

Returns

None

Return type

N/A

MOVE_CheckPlateWithTwoChannels(*device ML_STAR, variable i_int_FrontMostChannel, sequence i_seq_PlateToCheck, variable i_flt_TapWidth*)

This function will take the front most channel and the channel behind it to tap the labware at the sequence position. This function requires both channels to either have a tip or tool loaded on them before calling.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_FrontMostChannel

The front-most channel being used to check plates

Params i_seq_PlateToCheck

The sequence position to perform a tap to check for plate existence

Params i_fit_TapWidth

The distance in mm for the channels to be separated before tapping

Returns

Boolean determining whether plate was found (hslTrue) or not (hslFalse)

Return type

Boolean

MOVE_InitDispenseDrive(*device ML_STAR, variable i_int_ChannelNumber*)

This function moves the dispense drive for the specified 1mL channel to its home position

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_ChannelNumber

The channel to initialise the dispense drive

Returns

None

Return type

N/A

MOVE_InitDispenseDrive_5mL(*device ML_STAR, variable i_int_ChannelNumber*)

This function moves the dispense drive for the specified 5mL channel to its home position

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_ChannelNumber

The channel to initialise the dispense drive

Returns

None

Return type

N/A

PLATE_STACK_CountPlateStacks(*device ML_STAR, sequence i_seq_PlateStack_Full, sequence o_seq_PlateStack_Count, variable o_int_PlateCount*)

This function will use the channels to measure the number of plates in a plate stack

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_seq_PlateStack_Full

The full sequence of the plate stack to measure

Params o_seq_PlateStack_Count

The sequence of the plate stack measured

Params o_int_PlateCount

The total number of plates measured in the plate stack

Returns

None

Return type

N/A

QUERY_GetChannelPosition(*device ML_STAR, variable i_int_ChNumber, variable o_flt_XCoord, variable o_flt_YCoord, variable o_flt_ZCoord*)

This function will return the current coordinate of the specified 1mL channel

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_ChNumber

The number of the channel whose position is being checked

Params o_flt_XCoord

The X coordinate of the channel

Params o_flt_YCoord

The Y coordinate of the channel

Params o_flt_ZCoord

The Z coordinate of the channel

Returns

None

Return type

N/A

QUERY_GetChannelPosition_5mL(*device ML_STAR, variable i_int_ChNumber, variable o_flt_XCoord, variable o_flt_YCoord, variable o_flt_ZCoord*)

This function will return the current coordinate of the specified 5mL channel

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_ChNumber

The number of the channel whose position is being checked

Params o_flt_XCoord

The X coordinate of the channel

Params o_flt_YCoord

The Y coordinate of the channel

Params o_flt_ZCoord

The Z coordinate of the channel

Returns

None

Return type

N/A

QUERY_GetTipPresentState(*device ML_STAR, variable i_int_ChNumber, variable o_bln_TipPresent*)

This function outputs true if a tip is loaded and false if a tip is not loaded on the specified 1mL channel

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_ChNumber

The number of the channel which is being checked

Params o_bln_TipPresent

A boolean output of whether the tip is present (hslTrue) or not (hslFalse)

Returns

None

Return type

N/A

QUERY_GetTipPresentState_5mL(*device ML_STAR, variable i_int_ChNumber, variable o_bln_TipPresent*)

This function outputs true if a tip is loaded and false if a tip is not loaded on the specified 5mL channel

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_ChNumber

The number of the channel which is being checked

Params o_bln_TipPresent

A boolean output of whether the tip is present (hslTrue) or not (hslFalse)

Returns

None

Return type

N/A

QUERY_GetTipVolume(*device ML_STAR, variable i_int_ChNumber, variable o_ft_MaxVolume, variable o_ft_CurrentChannelVolume*)

This function queries the specified 1mL channel to get the max and the current channel volume. This volume includes the air gap and the conversion made by the correction curve.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_ChNumber

The number of the channel which is being checked

Params o_ft_MaxVolume

The maximum volume of the channel

Params o_ft_CurrentChannelVolume

The current volume in the channel

Returns

None

Return type

N/A

QUERY_GetTipVolume(*device ML_STAR, variable i_int_ChNumber, variable o_ft_MaxVolume, variable o_ft_CurrentChannelVolume*)

This function queries the specified 5mL channel to get the max and the current channel volume. This volume includes the air gap and the conversion made by the correction curve.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_int_ChNumber

The number of the channel which is being checked

Params o_ft_MaxVolume

The maximum volume of the channel

Params o_fit_CurrentChannelVolume

The current volume in the channel

Returns

None

Return type

N/A

SPLIT_WELLS_AddContainersToWell(*device ML_STAR, sequence i_seq_SequenceToSplit, variable i_int_SequenceIndex, variable i_int_MaxSplitNumber, sequence io_seq_SplitSequence*)

This function will split a well into a sequence of multiple containers, each of which can be aspirated from individually.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_seq_SequenceToSplit

The sequence that contains the well to be split into multiple containers

Params i_int_SequenceIndex

The index of the sequence that is to be split into multiple containers

Params i_int_MaxSplitNumber

The number of times that the selected container will be split. Won't exceed the width of the well

Params io_seq_SplitSequence

The sequence containing the split wells. Will append to the end of the sequence.

Returns

None

Return type

N/A

SPLIT_WELLS_Remove_Containers(*device ML_STAR, variable i_bool_UpdateVolumes*)

This function will remove the containers added by the split wells function

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_bool_UpdateVolume

A boolean determining whether to update the split source volume with the volume set in sample tracking.

Returns

None

Return type

N/A

TRAVEL_LANES_MoveChannelsToTravelLanes(*device ML_STAR*)

This function will move the 1mL channels to the travel lanes

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Returns

None

Return type

N/A

TRAVEL_LANES_MoveChannelsToTravellanes_5mL(*device ML_STAR*)

This function will move the 5mL channels to the travel lanes

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Returns

None

Return type

N/A

TRAVEL_LANES_MoveChannelsToYPosition(*device ML_STAR, sequence i_seq_TargetSequence, variable i_flt_XOffsetFromOrigin*)

Parameters include the Instrument, a destination sequence, and whether or not a shift in the x-direction is wanted before moving the 1mL channels in y-direction (to doubly ensure no crossver of open wells). Channels will be moved to their Y coordinates at the X origin of the next labware in the sequence plus the X offset. The channels will move to the current position of the input Sequence.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_seq_TargetSequence

The sequence to move the channels to

Params i_flt_XOffsetFromOrigin

The distance from the plate's origin to start the channel approach

Returns

None

Return type

N/A

TRAVEL_LANES_MoveChannelsToYPosition_5mL(*device ML_STAR, sequence i_seq_TargetSequence, variable i_flt_XOffsetFromOrigin*)

Parameters include the Instrument, a destination sequence, and whether or not a shift in the x-direction is wanted before moving the 5mL channels in y-direction (to doubly ensure no crossver of open wells). Channels will be moved to their Y coordinates at the X origin of the next labware in the sequence plus the X offset. The channels will move to the current position of the input Sequence.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_seq_TargetSequence

The sequence to move the channels to

Params i_flt_XOffsetFromOrigin

The distance from the plate's origin to start the channel approach

Returns

None

Return type

N/A

TRAVEL_LANES_MoveChannelsWithTravellanes(*device ML_STAR, sequence i_seq_TargetSequence, variable i_flt_XOffsetFromOrigin*)

This command is designed to shift all the 1mL channels on the instrument to either the front and/or the back of the instrument in a layout that ensures tips will not crossover any carriers. Parameters include the Instrument, a destination sequence, and whether or not a shift in the x-direction is wanted before moving the channels in

y-direction (to doubly ensure no crossover of open wells). Channels will be moved to their Y coordinates at the X origin of the next labware in the sequence plus the X offset. The channels will move to the current position of the input Sequence.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_seq_TargetSequence

The sequence to move the channels to

Params i_fit_XOffsetFromOrigin

The distance from the plate's origin to start the channel approach

Returns

None

Return type

N/A

TRAVEL_LANES_MoveChannelsWithTravelLanes_5mL(*device ML_STAR, sequence i_seq_TargetSequence, variable i_fit_XOffsetFromOrigin*)

This command is designed to shift all the 5mL channels on the instrument to either the front and/or the back of the instrument in a layout that ensures tips will not crossover any carriers. Parameters include the Instrument, a destination sequence, and whether or not a shift in the x-direction is wanted before moving the channels in y-direction (to doubly ensure no crossover of open wells). Channels will be moved to their Y coordinates at the X origin of the next labware in the sequence plus the X offset. The channels will move to the current position of the input Sequence.

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_seq_TargetSequence

The sequence to move the channels to

Params i_fit_XOffsetFromOrigin

The distance from the plate's origin to start the channel approach

Returns

None

Return type

N/A

TRAVEL_LANES_SingleSource_ChannelDisplacement(*device ML_STAR, variable i_strStepReturn*)

This command is designed to move unused 1mL channels to the back of the instrument when pipetting one well at a time while more than one channel has tips/liquid. It requires the Instrument type and Step Return variable from the Aspirate/Dispense step (this is used to determine which channel needs to be moved).

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_strStepReturn

The bound step return from given desired pipetting step

Returns

None

Return type

N/A

TRAVEL_LANES_SingleSource_ChannelDisplacement_5mL (*device ML_STAR, variable i_strStepReturn*)

This command is designed to move unused 5mL channels to the back of the instrument when pipetting one well at a time while more than one channel has tips/liquid. It requires the Instrument type and Step Return variable from the Aspirate/Dispense step (this is used to determine which channel needs to be moved).

Params ML_STAR

The ML_STAR itself, will be the only option in the dropdown

Params i_strStepReturn

The bound step return from given desired pipetting step

Returns

None

Return type

N/A

STRTOKENIZE

<https://github.com/theonetrue/nerd/VenusPackages/blob/main/StrTokenize.pkg>

The string tokenize library allows you to split a string into an array based on delimiters. It adds one function:

- *StrTokenize()*

StrTokenize(*variable strIn, variable strDelimiter, array arrTokens, boolean bAttendEmptyTokens*)

This function takes the input string and splits it at any point where it finds the specified delimiter, returning the substrings as members of an array. Can be told to add empty strings or not (occurs when the delimiter appears twice in a row).

Params strIn

The input string to be split

Params strDelimiter

The character to be used as the delimiter

Params arrTokens

The output array of split-substrings

Params bAttendEmptyTokens

Boolean determining whether to include empty substrings or not in the array

Returns

None

Return type

N/A

STRING (FROM HSLEXTENSIONS)

<https://github.com/theonetrue nerd/VenusPackages/blob/main/String.pkg>

The String library from HSLExtensions adds a few functions to facilitate easier manipulation of strings. It adds the following functions:

- `ConvertToAsciiArray()`
- `ConvertToCharArray()`
- `FromAsciiArray()`
- `Join()`
- `JoinWithDelimiter()`
- `Split()`
- `Trim()`

ConvertToAsciiArray(*variable i_strValue*)

This function converts the input string into an array with the regarding ASCII codes. If the input parameter is not a string the function returns an empty array.

Params i_strValue

The input string to be converted into an array

Returns

An array of ASCII codes, or an empty array if the input parameter is not a string

Return type

Array

ConvertToCharArray(*variable i_strValue*)

This function converts the input string into an array with the regarding characters. If the input parameter is not a string the function returns an empty array.

Params i_strValue

The input string to be converted into an array

Returns

An array of characters (strings with length 1), or an empty array if the input parameter is not a string

Return type

Array

FromAsciiArray(*array i_arrAsciiValues*)

This function converts an input array with ASCII codes into a string. If the input parameter is not an array with ASCII codes, the function returns an empty string.

Params *i_arrAsciiValues*

The input array of ASCII codes to be converted into a string

Returns

The output string formed by the concatenation of the converted versions of the ASCII codes. An empty string if the input parameter is not an array with ASCII codes.

Return type

Variable

Join(*array i_arrValues*)

This function joins an array of strings into a single string. Can be used to concatenate any number of strings into a single one. If the input parameter is not an array with strings, the function returns an empty string.

Params *i_arrValues*

The array of strings to be concatenated

Returns

The concatenated form of all the strings in the array, or an empty string if the input parameter is not an array of strings

Return type

Variable

JoinWithDelimiter(*array i_arrValues, variable i_strDelimiter*)

This function joins an array of strings into a single string and adds a delimiter between each substring. If the input parameter is not an array with strings, the function returns an empty string.

Params *i_arrValues*

The input array of strings to be concatenated

Params *i_strDelimiter*

The delimiter to be inserted between each substring

Returns

The concatenated strings from the array, with delimiters between each substring. An empty string if the input parameter is not an array of strings.

Return type

Variable

Split(*variable i_strValue, variable i_strDelimiter, variable i_bTrimWhitespaces*)

This function splits a string into substrings, forming an array of strings. The input string is split based on a delimiter that the user inputs.

Params *i_strValue*

The input string to be split into substrings

Params *i_strDelimiter*

The delimiter to be used to split the string

Params *i_bTrimWhitespaces*

Boolean determining whether leading and trailing whitespaces will be removed or not

Returns

An array of strings containing each substring formed from splitting the original string

Return type

Array

Trim(*variable i_strValue*)

This function trims leading and trailing whitespace characters from the input string. If the input parameter is not a string the function returns an empty string.

Params i_strValue

The input string to trim

Returns

The trimmed string

Return type

Variable

TOOLS LIBRARY

https://github.com/theonetrue nerd/VenusPackages/blob/main/Tools_Library.pkg

The tools library is designed to add some QoL functions to Venus. It adds the following functions.

- *CalculateNumberOfTubesAndVolumePerTube()*
- *Reformat_Sequence()*
- *Reformat_ValuesForDialog()*

CalculateNumberOfTubesAndVolumePerTube(*variable i_str_ReagentName, variable i_int_NumberOfSamples, variable i_ft_VolumePerSample, variable i_ft_MaximumContainerVolume, variable i_ft_DeadLabwareVolume, variable i_ft_DeadVolumeReagent, variable o_int_NumberOfContainersNeeded, array o_arr_VolumesPerContainer, array o_arr_DescriptionForDialog*)

This function calculates the number of containers needed for a reagent and the volume per container.

- params i_str_ReagentName**
The name of the reagent being calculated
- params i_int_NumberOfSamples**
The number of samples requiring the reagent
- params i_ft_VolumePerSample**
The volume of reagent required for each sample
- params i_ft_MaximumContainerVolume**
The maximum volume a single container can hold
- params i_ft_DeadVolumeLabware**
The dead volume in the labware
- params i_ft_DeadVolumeReagent**
The percentage dead volume of the reagent
- params o_int_NumberOfContainersNeeded**
The calculated number of containers needed
- params o_arr_VolumesPerContainer**
An array of the volumes required for each tube
- params o_arr_DescriptionForDialog**
An array of the descriptions for the dialog

Returns

None

Return type

N/A

Reformat_Sequence(*sequence io_sequenceInput, variable i_intNumberOfPositions*)

This function reformats a sequence to have as many positions as needed: Source Sequence: A1, B1

Total Number of positions needed = 4

Output Sequence: A1, B1, A1, B1

Params io_sequenceInput

The sequence of interest

Params i_intNumberOfPositions

The number of positions required

Returns

None

Return type

N/A

Reformat_ValuesForDialog(*variable i_flt_valueToReformat, variable o_str_valueFormatted*)

This function translates floats into strings and gets rid of any trailing 0s. Use the output value to display volumes in dialogs only.

Params i_flt_valueToReformat

The float to turn into a string

Params o_str_valueFormatted

The output string after it has been converted from a float

Returns

None

Return type

N/A

TRACELEVEL

<https://github.com/theonetrue nerd/VenusPackages/blob/main/TraceLevel.pkg>

The TraceLevel library adds a variety of functions related to the trace file. The functions it adds are:

- `GetTraceLevel()`
- `SetTraceLevel()`
- `SetStringIndicator()`
- `Trace_02()`
- `Trace_04()`
- `Trace_06()`
- `Trace_08()`
- `Trace_10()`
- `TraceArray()`
- `TraceArrayHorizontally()`
- `TraceArraysFaceToFace()`
- `TraceSequence()`
- `TraceSequenceParameter()`
- `TraceSequencePositions()`
- `TraceAction()`
- `SetActionIndicator()`

GetTraceLevel()

This function returns the current trace level.

Returns

This function returns one of the following predefined constants; `TRACE_LEVEL_NONE` (0), which corresponds to no traces at all. `TRACE_LEVEL_RELEASE` (1), which corresponds to only items with release trace level being traced. `TRACE_LEVEL_DEBUG` (2), which corresponds to everything being traced.

Return type

Variable

SetTraceLevel(*variable i_intTraceLevel*)

This function is used to set the trace level for the method or library.

Params i_intTraceLevel

The trace level for the library. Set to one of the following constants; TRACE_LEVEL_NONE (0), which corresponds to no traces at all. TRACE_LEVEL_RELEASE (1), which corresponds to only items with release trace level being traced. TRACE_LEVEL_DEBUG (2), which corresponds to everything being traced.

Returns

None

Return type

N/A

SetStringIndicator(*variable i_strStringIndicatorCharacter*)

This function is used to set one or more characters to indicate strings (by flanking them) in all traces. It can be useful to identify leading or trailing spaces in strings. You can do this by setting i_strStringIndicatorCharacter to ' or *.

Params i_strStringIndicatorCharacter

Characters that are being added at the beginning and the end of each string trace.

Returns

None

Return type

N/A

Trace_02(*variable i_intTraceLevel, variable i_varToTrace_01, variable i_varToTrace_02*)

This function is used to trace the value of 2 variables in one line. Will not automatically insert a space between the two variables.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_varToTrace_01

The first variable to trace.

Params i_varToTrace_02

The second variable to trace.

Returns

None

Return type

N/A

Trace_04(*variable i_intTraceLevel, variable i_varToTrace_01, variable i_varToTrace_02, variable i_varToTrace_03, variable i_varToTrace_04*)

This function is used to trace the value of 4 variables in one line. Will not automatically insert a space between the variables.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_varToTrace_01

The first variable to trace.

Params i_varToTrace_02

The second variable to trace.

Params i_varToTrace_03

The third variable to trace.

Params i_varToTrace_04

The fourth variable to trace.

Returns

None

Return type

N/A

Trace_06(*variable i_intTraceLevel, variable i_varToTrace_01, variable i_varToTrace_02, variable i_varToTrace_03, variable i_varToTrace_04, variable i_varToTrace_05, variable i_varToTrace_06*)

This function is used to trace the value of 6 variables in one line. Will not automatically insert a space between the variables.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_varToTrace_01

The first variable to trace.

Params i_varToTrace_02

The second variable to trace.

Params i_varToTrace_03

The third variable to trace.

Params i_varToTrace_04

The fourth variable to trace.

Params i_varToTrace_05

The fifth variable to trace.

Params i_varToTrace_06

The sixth variable to trace.

Returns

None

Return type

N/A

Trace_08(*variable i_intTraceLevel, variable i_varToTrace_01, variable i_varToTrace_02, variable i_varToTrace_03, variable i_varToTrace_04, variable i_varToTrace_05, variable i_varToTrace_06, variable i_varToTrace_07, variable i_varToTrace_08*)

This function is used to trace the value of 8 variables in one line. Will not automatically insert a space between the variables.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_varToTrace_01

The first variable to trace.

Params i_varToTrace_02

The second variable to trace.

Params i_varToTrace_03

The third variable to trace.

Params i_varToTrace_04

The fourth variable to trace.

Params i_varToTrace_05

The fifth variable to trace.

Params i_varToTrace_06

The sixth variable to trace.

Params i_varToTrace_07

The seventh variable to trace.

Params i_varToTrace_08

The eighth variable to trace.

Returns

None

Return type

N/A

Trace_10(*variable i_intTraceLevel, variable i_varToTrace_01, variable i_varToTrace_02, variable i_varToTrace_03, variable i_varToTrace_04, variable i_varToTrace_05, variable i_varToTrace_06, variable i_varToTrace_07, variable i_varToTrace_08, variable i_varToTrace_09, variable i_varToTrace_10*)

This function is used to trace the value of 10 variables in one line. Will not automatically insert a space between the variables.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_varToTrace_01

The first variable to trace.

Params i_varToTrace_02

The second variable to trace.

Params i_varToTrace_03

The third variable to trace.

Params i_varToTrace_04

The fourth variable to trace.

Params i_varToTrace_05

The fifth variable to trace.

Params i_varToTrace_06

The sixth variable to trace.

Params i_varToTrace_07

The seventh variable to trace.

Params i_varToTrace_08

The eighth variable to trace.

Params i_varToTrace_09

The ninth variable to trace.

Params i_varToTrace_10

The tenth variable to trace.

Returns

None

Return type

N/A

TraceArray(*variable i_intTraceLevel, variable i_strDescription, array i_arrvarToTrace*)

This function is used to trace an array of variables. It will trace each value of the array in its own line, along with the array description and the index of the value.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_strDescription

A description of the array, which will be at the start of each line of the array trace.

Params i_arrvarToTrace

The array to be traced.

Returns

None

Return type

N/A

TraceArrayHorizontally(*variable i_intTraceLevel, variable i_strDescription, array i_arrvarToTrace*)

This function is used to trace an array of variables. It will trace the array description, followed by each array index and value pair, all on one line.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_strDescription

A description of the array, which will be at the start of the array trace.

Params i_arrvarToTrace

The array to be traced.

Returns

None

Return type

N/A

TraceArraysFaceToFace(*variable i_intTraceLevel, variable i_strDescription_1, variable i_strDescription_2, array i_arrvarToTrace_1, array i_arrvarToTrace_2*)

This function is used to trace two arrays of variables at the same time, with values at the same index being shown on the same line as one another.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_strDescription_1

A description of the first array, which will be at the start of the array trace.

Params i_strDescription_2

A description of the second array, which will be at the start of the array trace.

Params i_arrvarToTrace_1

The first array to be traced.

Params i_arrvarToTrace_2

The second array to be traced.

Returns

None

Return type

N/A

TraceSequence(*variable i_intTraceLevel, sequence i_seqToTrace*)

This function is used to trace a sequence. It will list the sequence name, current position, the count and total positions in the sequence, the max number of positions available, and the number of used positions. It will then list the labware ID and position ID for each value of the sequence.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_seqToTrace

The sequence to be traced

Returns

None

Return type

N/A

TraceSequenceParameter(*variable i_intTraceLevel, sequence i_seqToTrace*)

This function is used to trace the parameters of a sequence. It will list the sequence name, current position, the count and total positions in the sequence, the max number of positions available, and the number of used positions.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_seqToTrace

The sequence to be traced

Returns

None

Return type

N/A

TraceSequencePositions(*device ML_STAR, variable i_intTraceLevel, sequence i_seqToTrace, variable i_blnCurrentPositionOnly*)

This function is used to trace the deck positions of a sequence. It will trace the sequence name, then for each part of the sequence it will trace the position ID and the x, y, z, and r coordinates of the position.

Params ML_STAR

The STAR device.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_seqToTrace

The sequence to be traced

Params i_blnCurrentPositionOnly

A boolean determining whether the function will trace positions for all sequence positions (0) or just the current sequence position (1).

Returns

None

Return type

N/A

TraceAction(*variable i_intTraceLevel, variable i_intAction, variable i_strFunctionName, variable i_strMethodName, variable i_strComment*)

This function is used to trace different action states for a module. [SetActionIndicator\(\)](#) can be used to help identify actions in log files more easily.

Params i_intTraceLevel

The trace level for the entry. Can be set to either TRACE_LEVEL_RELEASE (1) or TRACE_LEVEL_DEBUG (2). If set to 1, the function will show up in the trace when the trace level is set to either TRACE_LEVEL_RELEASE or TRACE_LEVEL_DEBUG, if set to 2, the function will only show up when the trace level is TRACE_LEVEL_DEBUG.

Params i_intAction

The action to be traced. Set to one of the following constants. START (1) corresponds to the entry being traced as starting. COMPLETE (2) corresponds to the entry being traced as finished successfully. ERROR (3) corresponds to the entry being traced as error occurred. PROGRESS (4) corresponds to the entry being traced as progressing. COMPLETE_WITH_ERROR (5) corresponds to the entry being traced as finished unsuccessfully.

Params i_strFunctionName

The function name for the action trace. Can be set to the return value of the HSL function [GetFunctionName](#), which will automatically be formatted correctly.

Params i_strMethodName

The method name for the action trace. Can be set to the return value of the HSL function [GetMethodFileName](#), which will automatically be formatted correctly.

Params *i_strComment*

A comment to be traced with the action trace.

Returns

None

Return type

N/A

SetActionIndicator(*variable i_intAction, variable i_strIndicator*)

This function is used to set one or more characters to indicate actions. By setting *i_strIndicator* to different characters you can make it easy to identify different actions in the trace.

Params *i_intAction*

The action to be traced. Set to one of the following constants. START (1) corresponds to the entry being traced as starting. COMPLETE (2) corresponds to the entry being traced as finished successfully. ERROR (3) corresponds to the entry being traced as error occurred. PROGRESS (4) corresponds to the entry being traced as progressing. COMPLETE_WITH_ERROR (5) corresponds to the entry being traced as finished unsuccessfully.

Params *i_strIndicator*

Character(s) that are being repeated at a total length of up to 100 and traced before and after the action trace.

Type

i_intAction: Variable

Returns

None

Return type

N/A

WINDOWS (FROM HSLEXTENSIONS)

<https://github.com/theonetrue nerd/VenusPackages/blob/main/Windows.pkg>

The windows library from HSLExtensions adds functions related to registries and special directories. The functions it adds are as follows:

- *GetRegistryValue()*
- *GetSpecialDirectory()*
- *SetRegistryValue()*

GetRegistryValue(*variable i_intHKey, variable i_strKey, variable i_strValueName, variable o_varValue*)

This function reads a value from the specified registry

Params i_intHKey

The main key of the registry. 1 = ClassesRoot, 2 = CurrentUser, 3 = LocalMachine, 4 = Users, 5 = CurrentConfig.

Params i_strKey

The key to the path (e.g. SOFTWAREPhoenixDirectoriesMethods)

Params i_strValueName

The name of the value to be read

Params o_varValue

The value which was read (type according to value type)

Returns

Boolean showing whether the key exists (hslTrue) or not (hslFalse)

Return type

Boolean

GetSpecialDirectory(*variable i_intSpecialDirectory*)

This function gets the specified special directory

Params i_intSpecialDirectory

The type of the special directory. 1 = WindowsDirectory, 2 = SystemDirectory, 3 = TemporaryDirectory.

Returns

The special directory

Return type

String

SetRegistryValue(*variable i_intHKey, variable i_strKey, variable i_strValueName, variable i_varValue, variable i_intValueType*)

This function writes a value to the specified registry

Params i_intHKey

The main key of the registry. 1 = ClassesRoot, 2 = CurrentUser, 3 = LocalMachine, 4 = Users, 5 = CurrentConfig.

Params i_strKey

The key to the path (e.g. SOFTWAREPhoenixDirectoriesMethods)

Params i_strValueName

The name of the value to be read

Params i_varValue

The value which was will be written in

Params i_intValueType

The type of the value. 1 = String, 2 = Number, 3 = Binary, 4 = ExpandableString

Returns

Boolean showing whether the function was successful (hslTrue) or not (hslFalse)

Return type

Boolean

ZERO UL SCANNER

<https://github.com/theonetrue nerd/VenusPackages/blob/main/ZerouLScanner.pkg>

The zero uL scanner library adds a single function which is designed to be added at the end of a method. It scans the trace file for any steps in which 0uL was pipetted and returns the number of times that happened. This is to check for steps where (for example) it is pipetting a volume based on a variable, and there is a spelling mistake in the variable, or the variable isn't defined correctly. It will work in both simulation mode and normal mode.

- *CheckFor0s()*

CheckFor0s()

This function checks for any steps in which 0uL has been pipetted and returns the number of times that has happened.

Returns

Number of times 0uL has been pipetted

Return type

Variable

A

AA_Abstract()
 built-in function, 45
 AddItemToArchive()
 built-in function, 77
 Append()
 built-in function, 25
 ArraySeqVLookup()
 built-in function, 31
 ArrayVLookup()
 built-in function, 31
 Assign()
 built-in function, 35

B

built-in function
 AA_Abstract(), 45
 AddItemToArchive(), 77
 Append(), 25
 ArraySeqVLookup(), 31
 ArrayVLookup(), 31
 Assign(), 35
 CalcAliquot_v2_1(), 21
 CalcAliquote(), 23
 CalcChannelPattern(), 22
 CalculateNumberOfTubesAndVolumePerTube(),
 119
 CHAN_ACCESS_Sort1Sequence(), 96
 CHAN_ACCESS_Sort1Sequence1Array(), 96
 CHAN_ACCESS_Sort1Sequence2Arrays(), 97
 CHAN_ACCESS_Sort2Sequences(), 98
 CHAN_ACCESS_Sort2Sequences1Array(), 99
 CHAN_ACCESS_Sort2Sequences2Arrays(), 100
 CheckFor0s(), 131
 CompareArrays(), 26
 Concat(), 26
 ContainsDuplicates(), 26
 ContainsValue(), 26
 ConvertArrayOfNumericIntegersToString(),
 32
 ConvertArrayOfNumericStringsToIntegers(),
 32

ConvertFileToASCII(), 37
 ConvertToAsciiArray(), 115
 ConvertToBooleanArray(), 27
 ConvertToCharArray(), 115
 ConvertToFloatArray(), 27
 ConvertToIntArray(), 27
 ConvertToStringArray(), 27
 Copy(), 28
 CopyActiveTraceFile(), 43
 CopyPlatePattern96ToTipRack(), 79
 CopyPlatePatternToPlate(), 79
 CreateMixedChannelPattern(), 91
 CreateOrderedChannelPattern(), 91
 FilCopyFileEx(), 61
 FilDeleteFileEx(), 61
 FilEof(), 55
 FilFindFile(), 56
 FilFindNextFile(), 56
 FilFormatBarcodeFile(), 56
 FilFormatReportFileEx(), 61
 FilGetBinPath(), 56
 FilGetCommState(), 57
 FilGetConfigPath(), 57
 FilGetLabwarePath(), 57
 FilGetLibraryPath(), 57
 FilGetLogFilesPath(), 57
 FilGetMethodsPath(), 58
 FilGetSystemPath(), 58
 FilIsNull(), 58
 FilReadString(), 58
 FilRemoveFields(), 58
 FilSearchPath(), 59
 FilSetCommState(), 59
 FilSetCommTimeouts(), 59
 FilUpdateRecord(), 59
 FilWriteString(), 60
 FindValue(), 28
 FromAsciiArray(), 115
 FullErrorAnalysis(), 39
 GenerateGeneralErrorsReport(), 43
 GeneratePipettingErrorsReport(), 43
 Get_ContainerBaseOffset(), 83

`Get_ContainerBaseThickness()`, 84
`get_distinct_from_array()`, 33
`Get_Height()`, 84
`Get_NameAndFileName()`, 84
`Get_NumberOfColumns()`, 85
`Get_NumberOfRows()`, 85
`Get_RackBaseToCoverBase()`, 85
`Get_StackHeight()`, 86
`Get_XY_dimensions()`, 86
`Get_XYZ_deckPosition()`, 86
`GetCarrierIDandSiteID_FromLabID()`, 83
`GetCommTimeouts()`, 57
`GetNumberOfPositionsLeft()`, 79
`GetRegistryValue()`, 129
`GetSpecialDirectory()`, 129
`GetTraceLevel()`, 121
`GetVersion()`, 51
`if_and()`, 81
`if_or()`, 81
`InitializeAllValues()`, 28
`InitializeArchive()`, 77
`IsBooleanArray()`, 28
`IsEmpty()`, 29
`IsFloatArray()`, 29
`IsIntegerArray()`, 29
`IsStringArray()`, 29
`Join()`, 116
`JoinWithDelimiter()`, 116
`LIQUID_LEVEL_GetLiquidLevelHeight()`, 101
`LIQUID_LEVEL_MeasureLiquidMulti()`, 102
`LIQUID_LEVEL_MeasureLiquidSingle()`, 102
`LIQUID_LEVEL_ReturnVolumesFromLiquidLevel()`, 103
`Lookup()`, 32, 89
`Make_Hidden()`, 41
`Make_ReadOnly()`, 41
`mergeArrays()`, 33
`MOVE_ChannelsToSequencePosition()`, 104
`MOVE_ChannelsToSequencePosition_5mL()`, 104
`MOVE_CheckPlateWithTwoChannels()`, 104
`MOVE_InitDispenseDrive()`, 105
`MOVE_InitDispenseDrive_5mL()`, 105
`PackArchive()`, 78
`ParseChannelPattern()`, 92
`Plate384PositionNumberToPositionText()`, 92
`Plate96PositionNumberToPositionText()`, 92, 93
`PLATE_STACK_CountPlateStacks()`, 105
`QUERY_GetChannelPosition()`, 106
`QUERY_GetChannelPosition_5mL()`, 106
`QUERY_GetTipPresentState()`, 106
`QUERY_GetTipPresentState_5mL()`, 107
`QUERY_GetTipVolume()`, 107
`Reformat_Sequence()`, 120
`Reformat_ValuesForDialog()`, 120
`Remove_Hidden()`, 41
`Remove_ReadOnly()`, 41
`removeValueFromArray_basedOnIndex()`, 34
`SeqEasyEdit()`, 80
`SetActionIndicator()`, 128
`SetRegistryValue()`, 129
`SetStringIndicator()`, 122
`SetTraceLevel()`, 51, 121
`Sort()`, 29
`Sort3ArraysByNumericAscendingOrder()`, 33
`Split()`, 116
`SPLIT_WELLS_AddContainersToWell()`, 108
`SPLIT_WELLS_Remove_Containers()`, 108
`Stat_Average()`, 63
`Stat_CorrelationCoefficient()`, 63
`Stat_Intercept()`, 64
`Stat_Min()`, 64
`Stat_RSQ()`, 64
`Stat_Slope()`, 64
`Stat_StdDeviation()`, 63
`STEP1_PrepareRegistryAndCfgFile()`, 45
`STEP2a_SimulateError_Channels()`, 45
`STEP2b_SimulateError_COREGripper()`, 46
`STEP2c_SimulateError_iSWAP()`, 46
`STEP2d_SimulateError_96Head()`, 47
`STEP2e_SimulateError_384head()`, 47
`STEP2f_SimulateError_BarcodeReading()`, 47
`STEP2g_SimulateError_Autoload()`, 48
`STEP2h_SimulateError_CRWashstation()`, 48
`STEP3_Restore_BackupCfgFile()`, 48
`STEP4_Optional_SwitchChecksum_ON()`, 49
`StrAsciiToStr()`, 68
`StrConcat12()`, 69
`StrConcat2()`, 68
`StrConcat4()`, 68
`StrConcat8()`, 68
`StrEvaluateExpr()`, 69
`StrFillLeft()`, 70
`StrFillRight()`, 70
`StrFind()`, 70
`StrFStr()`, 70
`StrFStrEx()`, 71
`StrFVal()`, 71
`StrGetLength()`, 71
`StrGetType()`, 71
`StrHexIStr()`, 71
`StrIsDigit()`, 72
`StrIsStr()`, 72
`StrIVal()`, 72
`StrLeft()`, 72
`StrMakeLower()`, 73

StrMakeLowerCopy(), 73
 StrMakeUpper(), 73
 StrMakeUpperCopy(), 73
 StrMid(), 73
 StrReplace(), 74
 StrReverseFind(), 74
 StrRight(), 74
 StrSpanExcluding(), 75
 StrStrToAscii(), 75
 StrTokenize(), 113
 StrTrimLeft(), 75
 StrTrimRight(), 75
 Trace_02(), 122
 Trace_04(), 122
 Trace_06(), 123
 Trace_08(), 123
 Trace_10(), 124
 TraceAction(), 127
 TraceArray(), 125
 TraceArrayHorizontally(), 125
 TraceArraysFaceToFace(), 125
 TraceSequence(), 126
 TraceSequenceParameter(), 126
 TraceSequencePositions(), 127
 TRAVEL_LANES_MoveChannelsToTravelLanes(), 108
 TRAVEL_LANES_MoveChannelsToTravelLanes_5mL(), 108
 TRAVEL_LANES_MoveChannelsToYPosition(), 109
 TRAVEL_LANES_MoveChannelsToYPosition_5mL(), 109
 TRAVEL_LANES_MoveChannelsWithTravelLanes(), 109
 TRAVEL_LANES_MoveChannelsWithTravelLanes_5mL(), 110
 TRAVEL_LANES_SingleSource_ChannelDisplacement(), 110
 TRAVEL_LANES_SingleSource_ChannelDisplacement_5mL(), 110
 Trim(), 116
 UnpackArchive(), 77
 Update_Value_in_Array(), 33
 UpdateLoadedLabware(), 54
 UpdateUsedLabware(), 53
 UpdateUsedPositions(), 53
 CalculateNumberOfTubesAndVolumePerTube()
 built-in function, 119
 CHAN_ACCESS_Sort1Sequence()
 built-in function, 96
 CHAN_ACCESS_Sort1Sequence1Array()
 built-in function, 96
 CHAN_ACCESS_Sort1Sequence2Arrays()
 built-in function, 97
 CHAN_ACCESS_Sort2Sequences()
 built-in function, 98
 CHAN_ACCESS_Sort2Sequences1Array()
 built-in function, 99
 CHAN_ACCESS_Sort2Sequences2Arrays()
 built-in function, 100
 CheckFor0s()
 built-in function, 131
 CompareArrays()
 built-in function, 26
 Concat()
 built-in function, 26
 ContainsDuplicates()
 built-in function, 26
 ContainsValue()
 built-in function, 26
 ConvertArrayOfNumericIntegersToString()
 built-in function, 32
 ConvertArrayOfNumericStringsToIntegers()
 built-in function, 32
 ConvertFileToASCII()
 built-in function, 37
 ConvertToAsciiArray()
 built-in function, 115
 ConvertToBooleanArray()
 built-in function, 27
 ConvertToCharArray()
 built-in function, 115
 ConvertToFloatArray()
 built-in function, 27
 ConvertToIntArray()
 built-in function, 27
 ConvertToStringArray()
 built-in function, 27
 Copy()
 built-in function, 28
 CopyActiveTraceFile()
 built-in function, 43
 CopyPlatePattern96ToTipRack()
 built-in function, 79
 CopyPlatePatternToPlate()
 built-in function, 79
 CreateMixedChannelPattern()
 built-in function, 91
 CreateOrderedChannelPattern()
 built-in function, 91

C

CalcAliquot_v2_1()
 built-in function, 21
 CalcAliquote()
 built-in function, 23
 CalcChannelPattern()
 built-in function, 22

F

`FilCopyFileEx()`
built-in function, 61

`FilDeleteFileEx()`
built-in function, 61

`FilEof()`
built-in function, 55

`FilFindFile()`
built-in function, 56

`FilFindNextFile()`
built-in function, 56

`FilFormatBarcodeFile()`
built-in function, 56

`FilFormatReportFileEx()`
built-in function, 61

`FilGetBinPath()`
built-in function, 56

`FilGetCommState()`
built-in function, 57

`FilGetConfigPath()`
built-in function, 57

`FilGetLabwarePath()`
built-in function, 57

`FilGetLibraryPath()`
built-in function, 57

`FilGetLogFilesPath()`
built-in function, 57

`FilGetMethodsPath()`
built-in function, 58

`FilGetSystemPath()`
built-in function, 58

`FilIsNull()`
built-in function, 58

`FilReadString()`
built-in function, 58

`FilRemoveFields()`
built-in function, 58

`FilSearchPath()`
built-in function, 59

`FilSetCommState()`
built-in function, 59

`FilSetCommTimeouts()`
built-in function, 59

`FilUpdateRecord()`
built-in function, 59

`FilWriteString()`
built-in function, 60

`FindValue()`
built-in function, 28

`FromAsciiArray()`
built-in function, 115

`FullErrorAnalysis()`
built-in function, 39

G

`GenerateGeneralErrorsReport()`
built-in function, 43

`GeneratePipettingErrorsReport()`
built-in function, 43

`Get_ContainerBaseOffset()`
built-in function, 83

`Get_ContainerBaseThickness()`
built-in function, 84

`get_distinct_from_array()`
built-in function, 33

`Get_Height()`
built-in function, 84

`Get_NameAndFileName()`
built-in function, 84

`Get_NumberOfColumns()`
built-in function, 85

`Get_NumberOfRows()`
built-in function, 85

`Get_RackBaseToCoverBase()`
built-in function, 85

`Get_StackHeight()`
built-in function, 86

`Get_XY_dimensions()`
built-in function, 86

`Get_XYZ_deckPosition()`
built-in function, 86

`GetCarrierIDandSiteID_FromLabID()`
built-in function, 83

`GetCommTimeouts()`
built-in function, 57

`GetNumberOfPositionsLeft()`
built-in function, 79

`GetRegistryValue()`
built-in function, 129

`GetSpecialDirectory()`
built-in function, 129

`GetTraceLevel()`
built-in function, 121

`GetVersion()`
built-in function, 51

I

`if_and()`
built-in function, 81

`if_or()`
built-in function, 81

`InitializeAllValues()`
built-in function, 28

`InitializeArchive()`
built-in function, 77

`IsBooleanArray()`
built-in function, 28

`IsEmpty()`

built-in function, 29
 IsFloatArray()
 built-in function, 29
 IsIntegerArray()
 built-in function, 29
 IsStringArray()
 built-in function, 29

J

Join()
 built-in function, 116
 JoinWithDelimiter()
 built-in function, 116

L

LIQUID_LEVEL_GetLiquidLevelHeight()
 built-in function, 101
 LIQUID_LEVEL_MeasureLiquidMulti()
 built-in function, 102
 LIQUID_LEVEL_MeasureLiquidSingle()
 built-in function, 102
 LIQUID_LEVEL_ReturnVolumesFromLiquidLevel()
 built-in function, 103
 Lookup()
 built-in function, 32, 89

M

Make_Hidden()
 built-in function, 41
 Make_ReadOnly()
 built-in function, 41
 mergeArrays()
 built-in function, 33
 MOVE_ChannelsToSequencePosition()
 built-in function, 104
 MOVE_ChannelsToSequencePosition_5mL()
 built-in function, 104
 MOVE_CheckPlateWithTwoChannels()
 built-in function, 104
 MOVE_InitDispenseDrive()
 built-in function, 105
 MOVE_InitDispenseDrive_5mL()
 built-in function, 105

P

PackArchive()
 built-in function, 78
 ParseChannelPattern()
 built-in function, 92
 Plate384PositionNumberToPositionText()
 built-in function, 92
 Plate96PositionNumberToPositionText()
 built-in function, 92, 93
 PLATE_STACK_CountPlateStacks()

built-in function, 105

Q

QUERY_GetChannelPosition()
 built-in function, 106
 QUERY_GetChannelPosition_5mL()
 built-in function, 106
 QUERY_GetTipPresentState()
 built-in function, 106
 QUERY_GetTipPresentState_5mL()
 built-in function, 107
 QUERY_GetTipVolume()
 built-in function, 107

R

Reformat_Sequence()
 built-in function, 120
 Reformat_ValuesForDialog()
 built-in function, 120
 Remove_Hidden()
 built-in function, 41
 Remove_ReadOnly()
 built-in function, 41
 removeValueFromArray_basedOnIndex()
 built-in function, 34

S

SeqEasyEdit()
 built-in function, 80
 SetActionIndicator()
 built-in function, 128
 SetRegistryValue()
 built-in function, 129
 SetStringIndicator()
 built-in function, 122
 SetTraceLevel()
 built-in function, 51, 121
 Sort()
 built-in function, 29
 Sort3ArraysByNumericAscendingOrder()
 built-in function, 33
 Split()
 built-in function, 116
 SPLIT_WELLS_AddContainersToWell()
 built-in function, 108
 SPLIT_WELLS_Remove_Containers()
 built-in function, 108
 Stat_Average()
 built-in function, 63
 Stat_CorrelationCoefficient()
 built-in function, 63
 Stat_Intercept()
 built-in function, 64
 Stat_Min()

built-in function, 64
 Stat_RSQ()
 built-in function, 64
 Stat_Slope()
 built-in function, 64
 Stat_StdDeviation()
 built-in function, 63
 STEP1_PrepereRegistryAndCfgFile()
 built-in function, 45
 STEP2a_SimulateError_Channels()
 built-in function, 45
 STEP2b_SimulateError_COREGripper()
 built-in function, 46
 STEP2c_SimulateError_iSWAP()
 built-in function, 46
 STEP2d_SimulateError_96Head()
 built-in function, 47
 STEP2e_SimulateError_384head()
 built-in function, 47
 STEP2f_SimulateError_BarcodeReading()
 built-in function, 47
 STEP2g_SimulateError_Autoload()
 built-in function, 48
 STEP2h_SimulateError_CRWashstation()
 built-in function, 48
 STEP3_Restore_BackupCfgFile()
 built-in function, 48
 STEP4_Optional_SwitchChecksum_ON()
 built-in function, 49
 StrAsciiToStr()
 built-in function, 68
 StrConcat12()
 built-in function, 69
 StrConcat2()
 built-in function, 68
 StrConcat4()
 built-in function, 68
 StrConcat8()
 built-in function, 68
 StrEvaluateExpr()
 built-in function, 69
 StrFillLeft()
 built-in function, 70
 StrFillRight()
 built-in function, 70
 StrFind()
 built-in function, 70
 StrFStr()
 built-in function, 70
 StrFStrEx()
 built-in function, 71
 StrFVal()
 built-in function, 71
 StrGetLength()

built-in function, 71
 StrGetType()
 built-in function, 71
 StrHexIStr()
 built-in function, 71
 StrIsDigit()
 built-in function, 72
 StrIsStr()
 built-in function, 72
 StrIVal()
 built-in function, 72
 StrLeft()
 built-in function, 72
 StrMakeLower()
 built-in function, 73
 StrMakeLowerCopy()
 built-in function, 73
 StrMakeUpper()
 built-in function, 73
 StrMakeUpperCopy()
 built-in function, 73
 StrMid()
 built-in function, 73
 StrReplace()
 built-in function, 74
 StrReverseFind()
 built-in function, 74
 StrRight()
 built-in function, 74
 StrSpanExcluding()
 built-in function, 75
 StrStrToAscii()
 built-in function, 75
 StrTokenize()
 built-in function, 113
 StrTrimLeft()
 built-in function, 75
 StrTrimRight()
 built-in function, 75

T

Trace_02()
 built-in function, 122
 Trace_04()
 built-in function, 122
 Trace_06()
 built-in function, 123
 Trace_08()
 built-in function, 123
 Trace_10()
 built-in function, 124
 TraceAction()
 built-in function, 127
 TraceArray()

- built-in function, [125](#)
- TraceArrayHorizontally()
 - built-in function, [125](#)
- TraceArraysFaceToFace()
 - built-in function, [125](#)
- TraceSequence()
 - built-in function, [126](#)
- TraceSequenceParameter()
 - built-in function, [126](#)
- TraceSequencePositions()
 - built-in function, [127](#)
- TRAVEL_LANES_MoveChannelsToTravelLanes()
 - built-in function, [108](#)
- TRAVEL_LANES_MoveChannelsToTravelLanes_5mL()
 - built-in function, [108](#)
- TRAVEL_LANES_MoveChannelsToYPosition()
 - built-in function, [109](#)
- TRAVEL_LANES_MoveChannelsToYPosition_5mL()
 - built-in function, [109](#)
- TRAVEL_LANES_MoveChannelsWithTravelLanes()
 - built-in function, [109](#)
- TRAVEL_LANES_MoveChannelsWithTravelLanes_5mL()
 - built-in function, [110](#)
- TRAVEL_LANES_SingleSource_ChannelDisplacement()
 - built-in function, [110](#)
- TRAVEL_LANES_SingleSource_ChannelDisplacement_5mL()
 - built-in function, [110](#)
- Trim()
 - built-in function, [116](#)

U

- UnpackArchive()
 - built-in function, [77](#)
- Update_Value_in_Array()
 - built-in function, [33](#)
- UpdateLoadedLabware()
 - built-in function, [54](#)
- UpdateUsedLabware()
 - built-in function, [53](#)
- UpdateUsedPositions()
 - built-in function, [53](#)